



**MALMET (AUSTRALIA)
PTY LTD**

Head Office and Customer Service
ABN 95 001 717 791

9-11 McKay Avenue
PO Box 373
Leeton NSW 2705

Phone: +61 2 6953 7677
Email: info@malmet.com.au

Bedpan / Urinal Bottle Washer Disinfector

Models ES915, ES935

Detergent Models with Hands Free Operation



Operation, Maintenance and Installation Manual

**Note: Due to Malmet's Policy of continuous product improvement;
design and technical specifications are subject to change without notice.**

| | |
|-----------------|---------------|
| Serial Number: | Supplied to: |
| Date Installed: | Installed by: |

It is important that the name from whom you purchased your device and the name of the installer are recorded above. The installer is responsible for the correct installation, start up and demonstrating the operation of this device. They are also responsible for issuing relevant certificates of compliance (these may differ from state to state).



Operation, Maintenance and Installation Manual

Table of Contents

| | |
|---|-----------|
| Foreword | 2 |
| Certifications and Compliances | 2 |
| Quality Policy | 2 |
| Important Warranty Reminder | 2 |
| Malmet Head Office and Factory Contact Details | 2 |
| Safety Instructions - Warnings | 3 |
| Intended Use | 4 |
| 1.0 Design Parameters | 4 |
| 1.1 Operating Cycle | 5 |
| 1.2 Chemical Dosing System | 6 |
| 1.3 Detergent | 6 |
| 1.4 Device Features | 8 |
| 1.5 Control Display Features | 8 |
| 1.6 Operating Features | 9 |
| 1.7 Device Labelling | 9 |
| 2.0 Installation and Commissioning | 10 |
| 2.1 Installation | 10 |
| 2.2 Service Connections | 12 |
| 2.3 Plumbing Connections | 12 |
| 2.4 Venting | 13 |
| 2.5 Electrical | 14 |
| 2.6 Service Connection Point | 15 |
| 2.7 Device Dimensions | 16 |
| 2.8 Installation clearances | 17 |
| 2.9 Commissioning (To only be completed by qualified persons) | 18 |
| 3.0 Loading | 19 |
| 4.0 Cycle of Operation | 21 |
| 4.1 Display Board | 26 |
| 4.2 Release of Processed items | 26 |
| 5.0 Maintenance | 27 |
| 5.1 Daily Maintenance (Operator or Maintenance Technician) | 27 |
| 5.2 Bi-Monthly Maintenance (Maintenance Technician) | 27 |
| 5.3 Recommended Preventative Maintenance Schedule | 28 |
| 5.4 Post Maintenance or Repair Safety Checks | 29 |
| 5.5 Door Safety Test | 29 |
| 5.6 Cycle Operation Check | 29 |
| 5.7 Post Maintenance and Repair of Electrical Systems | 29 |
| 5.8 Replacement of Safety Devices | 29 |
| 5.9 Validation | 29 |
| 5.10 Requalification | 30 |
| 5.11 Device faults | 30 |
| 5.12 Fault Indication | 30 |
| 5.13 Fault List | 31 |
| 5.14 Cold Water Overflow Indication | 33 |
| 5.15 Additional Faults not detected by micro-processor | 33 |
| 5.16 Making equipment safe after incomplete operating cycle | 34 |
| 6.0 Technical Specifications | 35 |
| 6.1 Power and Water Consumption | 35 |
| 6.2 Device Specifications | 35 |
| 6.3 Electrical Diagram – 20 Amp Models (ES915 Manufactured Post 28/02/2025) | 37 |
| 6.4 Electrical Diagram – 10 Amp Models (ES935 Manufactured Post 28/02/2025) | 38 |
| 6.5 Electrical Diagram – 20 Amp Models (ES915 Manufactured Pre 28/02/2025) | 39 |
| 6.4 Electrical Diagram – 10 Amp Models (ES935 Manufactured Pre 28/02/2025) | 40 |
| 6.5 Hydraulic Circuit Diagram - ES915/935 | 41 |
| Warranty Statement | 42 |



Operation, Maintenance and Installation Manual

Foreword

In order to obtain maximum life and efficiency from your Malmet Bedpan / Urinal Bottle Washer Disinfector and to aid in the safe operation of the device, please read and understand this manual thoroughly and follow all instructions before operating.

This manual provides information on the operation of the device. It is recommended that all persons operating the device have access to this manual for training purposes.

This device is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge.

Children should be supervised to ensure that they do not play with the device.

The specifications supplied in this manual were in effect at the time of publication. However, owing to Malmet (Australia)'s policy of continuous improvement, changes to these specifications may be made at any time without notice on the part of Malmet (Australia).

Certifications and Compliances

ARTG Identifier: 232012

Electrical Safety: Cert No. CS10847N to IEC 60601-1-2:2014 + A1:2020 (ED4.1)

EMC Compliance: Cert No. S2307003-2 V1

Watermark: Cert No. WMKA21156

Quality Policy

Malmet's quality management system is certified to ISO 13485 and ISO 9001 and guarantees the quality of this product.

Important Warranty Reminder

Should you have any problems with your device, contact the company from whom you purchased it, or Malmet (Australia) Pty Ltd.

It is important that the name from whom you purchased your device and the name of the installer are recorded on the front page of this manual. The installer is responsible for the correct installation, start up and demonstrating the operation of this device. They are also responsible for issuing relevant certificates of compliance (these may differ from state to state).

Malmet Head Office and Factory Contact Details

Malmet (Australia) Pty Ltd

9-11 McKay Avenue
PO Box 373
LEETON NSW 2705

Telephone: +61 2 6953 7677

E-mail: info@malmet.com.au

Website: www.malmet.com.au



Operation, Maintenance and Installation Manual

Safety Instructions - Warnings

Please read and understand this manual before using this device, if this device is used in a manner not specified by the manufacturer protection by the device may be impaired.

Please refer to this manual for information wherever this warning symbol is displayed -



Be aware of 240V Voltage



Disconnect power when servicing



Mains power ISO switch must be in an accessible position easily reached and not obstructed so device can be isolated from mains power during service



For the safe use of this device the responsible body should ensure that all operators are adequately trained to operate and maintain the device in its safe use



Failure to maintain device as per the requirements of this manual may impair performance



Only process items listed as per the design parameters in this manual



Be aware of steam discharge



Utensils and racks are hot to handle



Safety gloves and goggles must be worn when changing detergent



Hands-free sensor emits Infrared Radiation, avoid looking directly into the beam



Safety clothing with reflective tape can activate the hands-free sensor when device is in standby mode



Be aware of hot surfaces, pipes and hoses from steam and hot water



Machine to be serviced and repaired by trained personnel



Safety critical devices must only be replaced by trained personnel



Plumbing service connection must comply with AS/NZS 3500.1 and AS/NZS 3500.2



Not suitable for use in the presence of flammable anaesthetic mixtures with air or nitrous oxide and mode of operation as continuous



Operation, Maintenance and Installation Manual

Intended Use

The Malmet Bedpan / Urinal Bottle Washer Disinfector is intended to be used by health service organisation personnel for the cleaning and disinfection of medical devices during an automatic cycle of the types intended to be re-used as listed in the design parameters section of this manual.

1.0 Design Parameters

The Malmet Bedpan / Urinal Bottle Washer Disinfector has been designed within the following parameters:

- a) To empty, clean and disinfect devices of the type listed below intended to be re-used and that have been soiled by human excreta. *
- b) A single bedpan with lid and two urinal bottles can be emptied, cleaned and disinfected during each automatic cycle.

The utensils that can be cleaned in the device are: -

- i) Standard size bed pans
- ii) Standard commode bowl
- iii) Standard male and female urinal bottles
- iv) Most plastic urinal bottles including male non-spill and female urinal bottle
- v) Small slipper pan
- vi) Large slipper pan
- c) The cradle is designed to ensure that utensils are not dislodged during the cleaning cycle; the contents are emptied during door closure.
- d) The chamber and door are self-cleaning and do not permit water or soil to remain after a properly completed cycle. Steam disinfection ensures all internal surfaces are totally clean and safe.
- e) The flush and clean stage:
 - i) Removes the soil
 - ii) Clears the trap
- f) A complete cycle is completed in approximately four (4) minutes for the 20 Amp Bedpan /Urinal Bottle Washer Disinfector and six (6) minutes and thirty (30) seconds for the 10 Amp Bedpan / Urinal Bottle Washer Disinfector

Note: Dependant on ambient air temperature.

Note: The first cycle will take longer as the device has to reach the operating water temperature of 90°C.

- g) Malmet Door Obstruction Feature

This is factory set to 12Nm \pm 10%; causing the door to re-open when an obstruction is detected. This is designed to protect the integrity of the device and the articles within.

*Non-flushable items such as Feminine Hygiene products, Incontinence or Absorbent pads & diapers, Condoms, Dental Floss, Paper Towel, Cigarette butts, medications or chemicals cannot be disposed of in Malmet Washer Disinfectors and should be disposed as per facility &/or local waste management requirements.



Operation, Maintenance and Installation Manual

1.1 Operating Cycle

Stage 1 Flushing / Cleaning

One rotating 180° back spray

Two fixed 60° top sprays

Two fixed urinal sprays

Rinse for 10 seconds

Detergent pump on for 7 seconds to add 30ml of concentrate to pump

Main pump on for 3 seconds to apply detergent

Pause for 40 seconds

Rinse for 10 seconds to wash detergent off

Stage 2 Thermal Disinfection

Steam is delivered into the sealed wash chamber

Steam continues to heat the surface temperature of utensils and chamber to a minimum of 90°C for 60 seconds

Stage 3 Cool Down Rinse

One 60° fine mist top spray delivers cold water direct from mains.

Utensils are cooled to 55°C for safe handling.

In accordance and complies with:

Cleaning Efficacy

AS 5369, ISO 15883-1, ISO/TS 15883-5

Thermal Disinfection

AS 5369, ISO 15883-1



Operation, Maintenance and Installation Manual

1.2 Chemical Dosing System

The chemical dosing system ensures the correct amount of process chemicals are delivered at the correct stage during each operating cycle. Chemical addition is controlled by the automatic cycle controller and will display a fault if the incorrect amount of chemical is injected. The dosing system is designed for use with Malmest Chemical; other chemicals may impair the operation of the system.

1.3 Detergent

Malmest Washer Disinfectors are designed and tested using Malmest specific detergent. To ensure intended performance and trouble-free operation use only Malmest branded detergent, failure to use correct detergent may void warranty. The detergent system delivers 30ml of concentrated detergent in the initial wash stage.

Changing the Detergent Bottle



WARNING!



SAFETY GLOVES AND GOGGLES MUST BE WORN WHEN CHANGING DETERGENT AND DISPOSING OF EMPTY CONTAINERS

Current Safety Data Sheet for Malmest detergent is available in .PDF format from Malmest's website. www.malmest.com.au

Only use Malmest approved detergent (See technical data for detergent details)

1. Pull handle on detergent door and open.
2. Unscrew cap and pull out with suction hose (let hose hang in detergent chamber).
3. Remove empty bottle and replace with full bottle.

Note: Leave cap on new bottle until in position.

4. Remove cap on new bottle and fit existing hose and cap,
make sure suction hose is at bottom of bottle
5. Close detergent door.
6. Restart device operation as normal.



Disposal of Empty Detergent Bottles

1. When handling empty containers treat as though container is full and wear appropriate PPE as per MSD requirements.
2. Triple rinse container and puncture base to render unusable.
3. Dispose of containers as per facility chemical waste handling procedures, and/or local government requirements.



Operation, Maintenance and Installation Manual

Detergent Out Condition - F 30 Fault

If insufficient detergent is available to deliver the required quantity, the device will display an F30 dEt fault. Access to the load will be restricted until the condition is rectified.

To rectify this fault, install a new detergent bottle as per the procedure in this manual “Changing the Detergent Bottle”.

On completion of changing the detergent bottle, press either the “Manual Door” or “Manual Start” button to initiate a priming sequence. If priming is successfully the fault will clear and door will open. Inspect the items, ensuring they have been cleaned sufficiently before unloading. It is recommended the items are reprocessed in the event of an F30 dEt Fault.

If priming is unsuccessful, access to the load will continue to be restricted and F 30 dEt fault displayed.



To avoid the detergent out condition, check the low-level mark on detergent door, replace detergent when the detergent level reaches this mark

Operation, Maintenance and Installation Manual

1.4 Device Features


For the safe use of this device the responsible body should ensure that all operators are adequately trained to operate and maintain the device in its safe use.

Hands Free Operation

Zero Contamination – Totally Hands Free

The Malmest Bedpan / Urinal Bottle Washer Disinfector is fully automatic.

Door open, door close and cycle start can be activated without touching the device.



Front View


Manual operation of the door and a manual start button are available.

To Operate




- To **Open Door** – **Break Sensor Beam**
- To **Close Door** – **Break Sensor Beam**
- To **Start Cycle** – **Break Sensor Beam once when green LED's flash within 8 seconds after door closure**

Note: If door is opened then closed and cycle has not commenced the beam must be broken twice.

1.5 Control Display Features




1.6 Operating Features

| | | |
|---|--------------|--|
|  | POWER | On/Off Standby |
|  | MANUAL START | Start Cycle/Close door and start cycle |
|  | MANUAL DOOR | To manually open and close door |

1.7 Device Labelling

- 91-1002 ES915S-D S Trap 1 ph 20 Amp



Made in Australia by
MALMET (AUSTRALIA) PTY LTD
9-11 McKay Ave Leeton NSW 2705
ABN 95 001 717 791
Head Office: +61 2 6953 7677

MODEL **ES915S-D**


BATCH #


PART # **91-1002**

CERT # **CS10847N**

PHASES **1** VOLTS **240** AMPS **20.0** HERTZ **50**

SERIAL #

 REFER TO MANUAL
FOR POTENTIAL
HAZARDS
10.02.2025



MODEL **ES915S-D**

PART # **91-1002**


SERIAL #

MODEL **ES915S-D**

PART # **91-1002**

SERIAL #

- 91-1022 ES935S-D S Trap 1 ph 10 Amp



Made in Australia by
MALMET (AUSTRALIA) PTY LTD
9-11 McKay Ave Leeton NSW 2705
ABN 95 001 717 791
Head Office: +61 2 6953 7677

MODEL **ES935S-D**


BATCH #


PART # **91-1022**

CERT # **CS10847N**

PHASES **1** VOLTS **240** AMPS **10.0** HERTZ **50**

SERIAL #

 REFER TO MANUAL
FOR POTENTIAL
HAZARDS
10.02.2025



MODEL **ES935S-D**

PART # **91-1022**

SERIAL #

MODEL **ES935S-D**

PART # **91-1022**

SERIAL #

Issue 8

Page 9

March 2025



Operation, Maintenance and Installation Manual

2.0 Installation and Commissioning

2.1 Installation

To avoid problems with this device, these Installation Guidelines should be followed.

Installations must be carried out by a qualified and licenced tradesperson.

The device must be affixed to the building structure to ensure mechanical stability prior to use. The device can be affixed to the building structure via the holes provided in the device feet, or use of an appropriately rated appliance safety strap.


Prior to installation of the device, services as noted are to be provided by the facility. It is not the responsibility of Malmet to provide these service connections.

⚠ Overcurrent protection device-

- A 10A circuit breaker or fuse for 10A rated models must be installed in the building installation.
- A 20A circuit breaker or fuse for 20A rated models must be installed in the building installation.

Risk Assessment

It is recommended a risk assessment is conducted by the user both prior to and after installation and any risks identified mitigated to an acceptable level using the hierarchy of control;



<https://commons.wikimedia.org/w/index.php?curid=55610678>

Handling



Weights of Device: Net: 92 kg Shipping: 110 kg Shipping with crating: 151 kg.

- Handling of the device to installation site must be with fork lift or hand pallet truck.
- Before unpacking device inspect carton for any damage relating to forklift forks and damage relating to device falling over or for evidence of top loading
- After unpacking the device, inspect all external panels for damage.
- Remove the 4 screws holding the device to the pallet.



Operation, Maintenance and Installation Manual

- Follow your internal manual handling guidelines to manoeuvre the device off the pallet. The device can then be placed into position by fork lift or hand pallet truck.



Disposal of Packaging

- Please dispose of packaging as per facility procedures or local government requirements.



Disposal of Medical Device

- Please dispose of medical device as per state environmental regulatory requirements.

2.1.1 Positioning the ES-D

| Model | Placement | Access Required | Device Dimensions | | |
|---|--------------|-----------------|-------------------|------------|------------|
| | | | Height (mm) | Width (mm) | Depth (mm) |
| ES915 models ES935 models | Freestanding | Both Sides | 1305 | 595 | 660 |
| Please allow sufficient room for servicing purposes. Recommended space requirements 200mm on either side and 150mm at the rear of the device. | | | | | |

New Buildings

- Service connections are normally pre-placed after planning and consultation with all interested parties. Installation is by connection to the services provided.
- As the soil line (sewerage outlet) is the least flexible of all the connections, this usually influences the decision as to where to place the Bedpan / Urinal Bottle Washer Disinfector. If an existing soil line can be utilised this will represent a cost saving.
- The Bedpan / Urinal Bottle Washer Disinfector is supplied with either a 'S' or 'P' Trap as nominated by the Purchaser. The 'S' Trap connects through the floor and the 'P' Trap connects through the back wall. The trap section is easily removed if the wrong trap has been ordered. Refer to Diagram C2 for trap connections.
- Potential electromagnetic or other interference between other EQUIPMENT and other devices can possibly affect the Infra-red hands-free operation sensor. It is advisable to check all the equipment and devices in the intended installation area that have infra-red operation. **Electromagnetic interference can be prevented by installing the device in non-patient areas of the facility (or similar).**





Operation, Maintenance and Installation Manual

2.2 Service Connections

| MODEL | COLD WATER | SOIL LINE | ELECTRICAL |
|-------|-----------------------|-----------------------|------------------------------------|
| ES915 | GB $\frac{3}{4}$ Male | 100mm 'S' or 'P' Trap | 240V 1 phase @ 20 Amps 50 hertz |
| ES935 | GB $\frac{3}{4}$ Male | 100mm 'S' or 'P' Trap | 240V 1 phase @ 10 Amps 50 hertz |

2.3 Plumbing Connections



These installation guidelines must be followed to ensure the device will operate as intended.

Installations must be carried out by a qualified, licensed tradesperson.

Water Supply – Note: Water Supply connections must comply with AS/NZS3500.1 and be Watermark Certified.

This device requires only a potable cold-water connection; the water entry point is on the rear right-hand side of the device when viewed from the front. Please refer to the Service connection layout diagram for exact location. The device can be connected to any potable mains pressure cold water supply as a back-flow prevention air gap is incorporated in the design of the water tank. Complies to AS 2845.2

Cold Water – Note: See device specifications table for water quality requirements

Pressure: 100-400kPa

Temperature: 15-25°C



Nominal Flow: 10L/Min

Connection: GB $\frac{3}{4}$ " BSP

- The water supply is to be connected to an isolating valve or cistern stopcock placed approximately 1200mm from the floor to the right-hand side of the device (preferably not behind the device).
- If inlet water flow pressure is higher than 400 KPa an inline pressure reducing valve should be fitted.

Note: Old hose sets should not be re-used; it is recommended new compliant hose sets are used on installation.

Waste/Soil Line Connection



Note: Bell Ends can be cut off to suit

Diagram C2

For reasonable connection working space allow 150-200mm from rear of device to wall.

The device is 600mm wide and the centre of the trap is 300mm from each side.



Operation, Maintenance and Installation Manual

FREE STANDING 'S' TRAP PIPE POSITIONING

The centre of the soil line to receive the 'S' Trap should be approximately 272mm from the back wall. To allow for normal recommended minimum service access, space soil line 600/700mm from side wall.

If space restrictions do not allow for recommended side service access, Malmet suggest preference be given to providing the most space available on the right-hand side as you look at the front of the device. This will ease any difficulty in servicing the steam tank element and probe.

FREE STANDING 'P' TRAP PIPE POSITIONING

The centre of the soil line to receive the 'P' Trap should be approximately 410mm from the floor when the device is positioned 150mm from the wall. Because this pipe is graded to 5° this measurement will vary as the device is installed closer or further away from the back wall.

- There is no high temperature water discharge from device so no special high temperature pipe work is necessary. Soil line must comply to AS/NZS 3500.2
- Soil line connection is by a pan collar or other preferred method. If the belled end on the polyethylene moulded trap is not required it can easily be cut off to provide a straight pipe connection.
- Level the device by using the flanged screw in legs and, if possible, maintain approx. 100mm floor clearance for ease of floor cleaning. Malmet recommends affixing some of the leg flanges to the floor via stainless self-tapping screws to prevent sideways movements and damage to services and soil line connections.
- The soil line should protrude from the floor or wall at a minimum of 100mm.

2.4 Venting

No external vent pipework is required as the device is designed to condensate all visible steam within the device.



Operation, Maintenance and Installation Manual

2.5 Electrical



These installation guidelines must be followed to ensure the device will operate as intended and must be carried out by a qualified, licensed electrical tradesperson. The device must be installed and serviced to national wiring rules AS/NZS 3000.

| | | | |
|----------------------------|-----------|------|-------|
| Model: ES915 models | 240V 50Hz | 1 Ph | 20Amp |
| ES935 models | 240V 50Hz | 1 Ph | 10Amp |

Note: Mains power connection



Overcurrent Protection device: A circuit breaker or fuse must be installed in the building installation for the device.

10 Amp Devices – To be hard wired to a lockout ISO switch or Circuit breaker.



Switch must be wired to 10A max circuit breaker at main switch board

20 Amp Devices – To be hard wired to a lockout ISO switch or Circuit breaker.



Switch must be wired to 20A circuit breaker at main switch board

- An ISO switch or circuit breaker must be included in the installation. (Not supplied by Malmet)
- Mains power ISO switch must be suitably located and easily reached, approximately 1500mm above floor level adjacent to device. The 1.7m mains power lead exits the device approximately 1200mm above floor level on the right-hand side of the device.
- It must be marked as the disconnecting device for the equipment.
- It must have contact disconnection of all poles to provide full disconnection.
- Ensure isolator is suitably placed so device will not obstruct safe access to the isolator once installed.
- If the supply lead is damaged, it shall be replaced by Malmet; its authorised service agent or similarly qualified person in order to avoid hazard.



Mains power ISO switch or circuit breaker must be in an accessible position so device can be isolated from mains power during service.


All devices must be earthed.



The maximum permissible system impedance Z_{sys} :
 $Z = 0.220 \text{ Ohm} + j 0.137 \text{ Ohm} (0.220 \text{ Ohm} + 437 \text{ } \mu\text{H})$

Operation, Maintenance and Installation Manual


2.6 Service Connection Point



Dimension Tolerance $\pm 5\text{mm}$

Operation, Maintenance and Installation Manual



2.7 Device Dimensions



Dimension Tolerance $\pm 5\text{mm}$

Operation, Maintenance and Installation Manual

2.8 Installation clearances



Dimension Tolerance $\pm 5\text{mm}$



Operation, Maintenance and Installation Manual

2.9 Commissioning (To only be completed by qualified persons)

- a) Before switching on the device make sure the **DEVICE IS LEVEL** and **WATER TAPS ARE ON**. Check that the **DRAIN WASTE** is connected.
- b) Turn on the power at the isolation switch and press the standby button on the front display. The digital display will illuminate.
- c) Check that the water tank is filling with cold water and that it fills to Level 3. Ensure the lid is put back on water tank.
- d) Check that the steam generator tank has filled to the full Level 3. The steam generator will preheat to 85°C in standby mode or during cycle.
- e) Check that the 5-litre detergent bottle has been fitted.
- f) Flush approximately ½ litre of water down the steam generator tank overflow pipe. This will fill the 'S' Trap at the hose junction and prevent steam coming back up into the water tank.
- g) Open the door and turn power off at the control. Check that flush nozzle rotate freely. Check that all spray nozzles are tight. Turn power back on and close door by pressing the manual door button.
- h) Re-install all covers
- i) The device will display 'rdY' (ready) when the water levels are reached.
- j) Run a cycle and conform cycle is complete.

Note: DO NOT USE THE DEVICE WITHOUT THE WATER SUPPLY TURNED ON

Operation, Maintenance and Installation Manual


3.0 Loading



For the safe use of this device the responsible body should ensure that all operators are adequately trained to operate and maintain the device in its safe use.

The device is intended to be loaded as per the below configurations, please refer to the design parameters section of this manual for items that this device is designed to process. Part loads or singular items can be processed but must be loaded in their designated positions as below.

Bed Pan + 2 Urinal Bottles



Slipper Pan + 2 Urinal Bottles



Push down slipper pan
until it locks into
position

Operation, Maintenance and Installation Manual

Standard Commode Bowl + 2 Urinal Bottles



Large Commode Bowl + 2 Urinal Bottles




Operation, Maintenance and Installation Manual

4.0 Cycle of Operation


1. Once the Bedpan / Urinal Bottle Washer Disinfector is installed into position with power and water connected, push power on button.
2. When the Bedpan / Urinal Bottle Washer Disinfector is turned on via the display panel the **display shows fill** indicating steam generator and cold-water tanks are filling to the high-water level.

At this stage the device will not commence a cycle until Step 5.



3. Once the steam tank water level probes sense the tank is full it will start heating.
4. The steam tank is heated to temperature 90°C - 92°C. The Bedpan / Urinal Bottle Washer Disinfector is now ready to commence a cycle.


The display shows rdY (ready).



5. Break the infra-red beam by moving your hand between the infra-red indicators, this will open the door.

The display shows open.


Note: Whenever the door is opened the elements are in the off mode.



Operation, Maintenance and Installation Manual

6. Once the Bedpan / Urinal Bottle Washer Disinfector is “loaded”, break the infra-red beam again to close the door.

The display shows Shut




7. When the door is fully closed, 2 green LED's will flash in the infra-red sensor area for 8 seconds. Break the beam again to start a cycle.

The display shows rdY

If the 8 seconds has elapsed the flashing LED's will stop and the Bedpan / Urinal Bottle Washer Disinfector will not start for a cycle. The operator has to repeat the process from Step 5 (or press manual start).

Note: If this occurs the beam must be broken twice.


Note: On 10A models, if the steam generator temperature is not at 90°C after the beam has been broken or the manual start button pressed, the display will show “CYCL”. The Green LED's will stop flashing and stay on. When the steam generator has reached temperature, the display will show rinS and the wash cycle will start automatically.



Operation, Maintenance and Installation Manual


8. Once a cycle has been initiated the elements are turned on in the steam tank to start heating the water.
The wash pump is turned on for 10 seconds for the initial wash.

The display shows rinS




9. Detergent pump comes on delivering 30ml of concentrate to rinse pump.

The display shows dEt




10. Rinse pump comes on to distribute detergent over and inside items

The display shows rinS



11. Pause for 40 seconds so detergent can take effect


The display shows PAUS




Operation, Maintenance and Installation Manual

12. Rinse pump comes on for 10 seconds to rinse detergent from items

The display shows rinS



13. Once the wash cycle is complete the display will begin to flash the internal chamber temperature during heating until the disinfection temperature is reached. At this point the fill solenoid is disabled and will not operate until Step 19.



14. The elements in the steam generator heat to produce steam, which is directed into the wash chamber raising the internal cabinet temperature above 90°C.
15. Once the internal chamber temperature has reached 90°C the temperature display will stop flashing and the disinfection cycle will commence.



16. All elements on 20 Amp models and 10 Amp models will continue to heat during the disinfection cycle until the internal chamber temperature reaches 91°C.
17. For 20 Amp models, at 91°C one element will turn off and at 92°C another element will turn off.
For 10 Amp models, at 91°C one element will turn off.

For all models, one element will remain on until 15 seconds before the end of the cycle. If the internal chamber temperature falls below 92°C, the second element will turn back on and if the chamber temperature falls below 91°C, the third element will turn back on.

Note: Regardless of element combinations 15 seconds before the end of the cycle ALL elements will turn off. This is to overcome the effects of overshoot in the steam generator producing steam at the end of the disinfection cycle.

18. The maximum time allowable for the disinfection cycle is 10 minutes. If this time is exceeded then the device will go to fault and elements turn off.




Operation, Maintenance and Installation Manual


19. After the disinfection cycle is complete the cooling solenoid will turn on for 15 seconds. This will cool down the surface temperature of the utensils and make them safe to handle.

The display shows cool.

This water is direct from the mains.



20. At the end of the cool down rinse, **the display shows End**. This will be displayed until the door is re-opened. If the chamber temperature is above the safe opening temperature the display shows "HOT" and access is restricted until temperature has reduced to a safe level.



21. If the Bedpan / Urinal Bottle Washer Disinfector is not used, the steam generator temperature will be kept at 85°C-90°C for a period of 15 minutes. After 15 minutes, if the Bedpan / Urinal Bottle Washer Disinfector is still not used it will enter into idle mode and the elements will only come on at 50°C and rise to the temperature of 90°C repeatedly until the cycle is started again.



Operation, Maintenance and Installation Manual

4.1 Display Board

| | | Green LED Display |
|---|-------------|----------------------|
| Door open | | OPEN |
| Door closed | | Shut |
| Steam generator fill | | Fill |
| Cycle ready | | rdY |
| Cycle started | | rinS |
| Detergent | | dEt |
| Pause | | PAUS |
| Disinfection cycle started | (flashing) | 90°C |
| Disinfection temperature reached | | 'actual' temp °C |
| Cool down final rinse | | cool |
| End of cycle | | End |
| Chamber above safe opening temperature | | Hot |
| Yellow LED illuminates on control panel | | Overflow |
| 2 flashing green LED in activation area | (8 seconds) | Cycle can be started |

4.2 Release of Processed items

On completion of the cycle, it is recommended that all processed items are unloaded and visually inspected as clean as per the requirements AS 5369. It is recommended all items are removed on completion of the cycle and stored in a designated area to prevent any risk of contamination or mixing with unprocessed items.



If the processing status of a load or item is unknown or in question, it should be treated as unprocessed and not released.



Operation, Maintenance and Installation Manual

5.0 Maintenance



All maintenance; preventative or breakdown shall be carried out by a qualified person. Failure to comply with this condition may result in unsafe conditions.

Preventative maintenance must be completed as per instructions below. Failure to comply may result in warranty claims being rejected.

The Malmet Bedpan / Urinal Bottle Washer Disinfector is self-cleaning, however proper care should be taken to ensure that the device is cleaned and maintained in accordance with maintenance instructions for Malmet Bedpan / Urinal Bottle Washer Disinfectors and in accordance with all other regulatory and common-sense practices.

Preventative Maintenance Schedule

5.1 Daily Maintenance (Operator or Maintenance Technician)

- Wipe out the inside of the door and chamber with warm water and detergent. A wipe with disinfectant is also desirable.
- Wipe over outside stainless steel panels with a stainless-steel cleaner.
- Wipe the front control panel with a soft cloth and mild detergent as necessary. Care should be taken not to damage the digital display or to activate a cycle.
- Check level in detergent bottle, replenish as necessary.
- Check Bedpan wash spinner turns freely by manually spinning nozzle.
- Check spray nozzles are not blocked, insert a small pin into the orifice of the spray nozzle to ensure no blockages are present. If a build-up or blockage is noticed please have a technician clean nozzle as per Bi-monthly maintenance to ensure adequate performance.
- Visually inspect for signs of leaking fluid from the Bedpan / Urinal Bottle Washer Disinfector, as these may pose a pathogenic risk. Inspect around the chamber door, beneath the device and service connection points. If a spill or leak is observed, clean up as per the facility's infection control procedure and notify a maintenance technician to fix the cause of the leak.

5.2 Bi-Monthly Maintenance (Maintenance Technician)



WARNING! 240 VOLTS

ISOLATE DEVICE FROM ELECTRICAL SUPPLY BEFORE SERVICING



HOT SURFACES!

ALLOW DEVICE TO COOL PRIOR TO COMMENCING SERVICE WORKS

- Inspect for steam or fluid leaks, tighten unions, hose clamps and glands where necessary.
- Check flush nozzles for free rotation.
- Check spray nozzles are not blocked, insert a small pin into the orifice of the spray nozzle to ensure no blockages are present. If a build-up or blockage is detected have technician clean nozzles as per nozzle cleaning instructions below to ensure adequate performance.
- Remove level probe in steam generator tank and clean off any build-up of residue.
- Remove temperature probes in steam tank and chamber and clean off any build-up of residue.
- Visually inspect build-up of residue in steam tank, especially in areas of poor water quality.
- Check filter in the water inlet solenoid valve and clean as necessary.
- Check all electrical connections, and tighten if necessary.



Operation, Maintenance and Installation Manual

Nozzle Cleaning Instructions

Remove and clean the sprays in sequence (do not mix sprays with other spray groups).

Chamber top to bottom: Cool down sprays x 1, Top sprays x 2, and rotary nozzles x 1,

Urine bottle sprays x 2 bottle hooks (unscrew from inside chamber)

Hold under tap and pressurise through the nozzle outlet in opposite direction of normal flow or clean in ultrasonic cleaner if available. An appropriately sized Oxy-Acetylene tip cleaner pin may be used to clean out spray nozzle orifices. In areas with hard water or high minerals in water supply, chemical de-scaler may be the quickest and easiest means to remove build-up of deposits. Replace before removing next spray group

Stainless Steel Maintenance/Care

Under normal usage, stainless steel products require regular cleaning with a soft clean rag moistened with a mild detergent followed by a water-moistened clean rag and then a dry rag.

The #4 satin finish steel should be protected against Muriatic acid and caustic or abrasive materials and harsh cleaning detergents. In the event such agents cause discoloration, polish with a stainless-steel cleaner such as 3M Stainless Steel Cleaner & Polish and a pad.

5.3 Recommended Preventative Maintenance Schedule

(to be performed by qualified maintenance personnel)

It is recommended that preventive maintenance is performed by a qualified maintenance technician every **3 Months** for devices in high use, **6 Months** for medium use or **Annually** for low use devices



WARNING! 240 VOLTS

ISOLATE DEVICE FROM ELECTRICAL SUPPLY BEFORE SERVICING



HOT SURFACES!

ALLOW DEVICE TO COOL PRIOR TO COMMENCING SERVICE WORKS

- a) Remove top and side panels and front bottom panel.

Note: Panel removal

- i) Remove 2 screws on the top panel.
 - ii) Push the side panel back then lift up to remove.
 - iii) Remove 2 screws from bottom front panel.
- b) Remove three sprays, two bottle hooks and rotary nozzle from inside the chamber – check that the holes are clear. Clean as required as per **5.2 Bi-monthly Maintenance Nozzle Cleaning Instructions**.
- c) Start cycle and check that cold water pump is working, 17 seconds of the cycle.
- d) After ascertaining that the pump is working properly, check the solenoid operation in the cold-water tank. Check steam sensor (steam generator) maintenance.
- e) Make sure the solenoids are completely shutting off and levels are not creeping up, if so, clean and/or replace the solenoid.
- f) Tighten the screws retaining the pan rack on the door.

Only if necessary after visual check

- g) Clean dust and grit off components.

Malmet will make available on request circuit diagrams, component parts lists, descriptions, calibration instructions, or information which will assist the user's appropriately qualified technical personnel to repair those parts of the device.



Operation, Maintenance and Installation Manual

5.4 Post Maintenance or Repair Safety Checks

After servicing or repair the device shall be checked to confirm correct cycle operation and that the device is in a safe operating state before being returned to service.

5.5 Door Safety Test

To check the correct functionality of the door safety and position sensors, open and close the door using the manual door open button. During operation the device will conduct a self-check for correct operation of these switches. Once complete, check door is firmly closed and is flush with outer panel and then press the cycle "Start" button. If the wash pump starts, the self-check has been completed successfully and the door safety devices are functioning correctly.

5.6 Cycle Operation Check

After completion of the required checks above, a cycle operation check is to be completed. Run the device through a complete cycle ensuring correct sequence of operation as per page 12 of this manual. Inspect for fluid leaks from around the chamber door, rear service connections and beneath the device. If no faults are found during the inspection, and the device displays no faults at end of the cycle the device is fit for return to operation.

5.7 Post Maintenance and Repair of Electrical Systems

Additional to the above safety tests if any repair or service work has been carried out that could affect the electrical safety of the device; it must be inspected and tested as per the requirements of AS/NZS3760 prior to being returned to service.

Information on replacement of these devices can be found in the Service Technicians Manual.

If device fails any of the tests it is not safe for operation, and cannot be put into service until the cause of the fault has been rectified and successful completion of tests.

5.8 Replacement of Safety Devices

The following safety devices must be replaced and tested by a Malmet trained technician, failure to do so may impair the protection by the device.

- Door Actuator
- Door Micro Switch
- Element Thermal Cut-out

Information on replacement of these devices can be found in the Service Technicians Manual.

5.9 Validation

It is recommended that Validation of the device is carried out in accordance with the requirements of AS 5369 & ISO 15883, consisting of 3 identified stages.

IQ – Installation Qualification

OQ – Operational Qualification

PQ – Performance Qualification

Installation Qualification

Performed upon installation of equipment to ensure equipment has been installed in accordance with applicable Standards and manufacturer's instructions.

Operation, Maintenance and Installation Manual

Operation Qualification

Performed immediately after installation, relocation, change of a service to the device or after repair prior to device being put back in service in accordance with applicable Standards.

Performance Qualification

Performed immediately after IQ or OQ for new installs or relocated equipment, or change of a service to the device or repair that might adversely impact the quality of the RMD.

5.10 Requalification

It is recommended requalification of the process is performed annually in accordance with the requirements of AS 5369 and ISO 15883 applicable parts

5.11 Device faults

In the event of a fault or malfunction with the device during operation, abort the cycle by pressing the Power button on the control LCD. Once the device is off, isolate the mains power supply and water services. Notify maintenance or service provider of the fault and ensure device is not used until the fault is rectified.

5.12 Fault Indication

The device is controlled by a micro-processor. The processor has fault detection capability and indicates faults by code on the digital display; Diagram B1 (below) shows an example of fault code displayed.




Diagram B1
(Fault 01 indicated)

Please see table in 5.13 Fault List for a full list of fault codes along with their description and possible cause.

Note: If a fault is displayed on the LED display, turn device off and turn it back on, press START to try to complete a cycle. If a fault persists contact a Service Technician. Record the fault code that has appeared for the Service Technician.

Table B1 indicates these faults. For a more detailed description of these fault conditions see Appendix A – Faults.



Operation, Maintenance and Installation Manual

5.13 Fault List

Following is list of fault conditions and most possible causes.

| | | |
|-----------------|---|--|
| Fault 01 | Door Open Fault – Not confirmed door closed (i.e. door switch still indicates open) | |
| | <i>Cause</i> | <ul style="list-style-type: none"> • Door limit switch stuck out / broken |
| | Macerator Pump Fault – Where Fitted (Macerator Pump high level switch activated) | |
| | <i>Cause</i> | <ul style="list-style-type: none"> • No Power to Macerator Pump • Pump outlet/waste blockage • Faulty or Stuck Level Sensor |
| Fault 02 | Wash / Disinfection Cycle Time Out (taking too long to complete cycle) | |
| | <i>Cause</i> | <ul style="list-style-type: none"> • Steam leak • Final rinse solenoid stuck open • Faulty temperature sensor - hot water tank • Blown element leg • Steam condensate leaking cold water into chamber • Faulty temperature sensor - chamber |
| Fault 03 | Faulty Hot Water Probes (high indicates water, but low does not indicate water) | |
| | <i>Cause</i> | <ul style="list-style-type: none"> • Faulty water level probe |
| Fault 04 | Door motor activation time out (taking too long to close or open door) | |
| | <i>Cause</i> | <ul style="list-style-type: none"> • Faulty door actuator • Door axle bar sticking |
| Fault 06 | Hot water tank fill time exceeded (high level not reached during pre-fill, resulting in time out) | |
| | <i>Cause</i> | <ul style="list-style-type: none"> • Faulty fill solenoid • Faulty water level probe • Cold water tank empty or not filling quick enough • Low water pressure |
| Fault 07 | Hot water tank low water level sensed during cycle | |
| | <i>Cause</i> | <ul style="list-style-type: none"> • Faulty water level probe • Cold water tank empty or not filling quick enough • Leaking tank or hose • Element stuck on • Condensing coil blocked |
| Fault 08 | Door closure fault – not confirmed door open (door switch still indicates closed) | |
| | <i>Cause</i> | <ul style="list-style-type: none"> • Door limit switch stuck on |
| Fault 09 | Failed to reach disinfection start temperature, resulting in time out | |
| | <i>Cause</i> | <ul style="list-style-type: none"> • Steam leak • Final rinse solenoid stuck open • Faulty temperature sensor hot water tank • Blown element leg • Steam condensate leaking cold water into chamber • Faulty temperature sensor chamber • Door not reaching interlocking micro switch • Door not fully closed • Over temperature board (power to element open circuit) • Element over temp protection manual reset thermal cut-out switch open circuit |



Operation, Maintenance and Installation Manual


| | | |
|------------------|---|--|
| Fault 10 | Disinfection temperature not maintained, resulting in cycle time out | |
| | Cause | <ul style="list-style-type: none"> • Steam leak • Final rinse solenoid stuck open • Faulty temperature sensor hot water tank • Blown element leg • Steam condensate leaking cold water into chamber • Faulty temperature sensor - chamber |
| Fault 12 | Hot tank short temperature sensor | |
| | Cause | <ul style="list-style-type: none"> • Faulty temperature sensor steam generator |
| Fault 13 | Hot tank open circuit temperature sensor | |
| | Cause | <ul style="list-style-type: none"> • Faulty temperature sensor steam generator |
| Fault 14 | Disinfection chamber short circuit temperature sensor | |
| | Cause | <ul style="list-style-type: none"> • Faulty temperature sensor - chamber |
| Fault 15 | Disinfection chamber open circuit temperature sensor | |
| | Cause | <ul style="list-style-type: none"> • Faulty temperature sensor - chamber |
| Fault 30 dEt | Insufficient detergent last wash | |
| | Cause | <ul style="list-style-type: none"> • Empty detergent bottle • Faulty detergent pump • Blocked detergent line |
| Fault 31 HEAt | Disinfection temperature not maintained on 2nd probe | |
| | Cause | <ul style="list-style-type: none"> • Steam leak • Final rinse solenoid stuck open • Faulty temperature sensor hot water tank • Blown element leg • Steam condensate leaking cold water into chamber • Faulty temperature sensor – chamber • Faulty temperature sensor – chamber 2 |
| Fault 32 CFSd | Unsupported fault (cycle did not follow supported order) | |
| | Cause | <ul style="list-style-type: none"> • Interrupted power • Processor glitch |
| Fault 33 bLOC | Chamber blockage present – flood sensors (Sensors are detecting a build-up above the bottom of the door inside the chamber) | |
| | Cause | <ul style="list-style-type: none"> • Blocked trap • Blocked / backed up drain • Dirty sensors |
| Fault 34 dEtH | Extra detergent injected last wash (the volume of detergent injected during the cycle was greater than necessary) | |
| | Cause | <ul style="list-style-type: none"> • Faulty flow meter • Faulty detergent pump relay • Failed check valve |
| Fault 35 PobE | Temperature discrepancy – temperature difference >5°C between probes | |
| | Cause | <ul style="list-style-type: none"> • Faulty temperature sensor/s • Fluctuating temperatures |
| Fault 38 FPri | Failed priming – exceeded allowed priming attempts | |
| | Cause | <ul style="list-style-type: none"> • Detergent bottle not changed • Detergent line split / leaking • Flow meter faulty • Check valve failed |

Operation, Maintenance and Installation Manual

5.14 Cold Water Overflow Indication

The yellow LED is used for cold water tank overflow.

Overflow is indicated by a STEADY ON LED, whenever the overflow level is detected.



Should the cold-water tank overflow into the soil line the water overflow YELLOW LED will illuminate.

5.15 Additional Faults not detected by micro-processor

These are faults that may occur which the micro-processor cannot detect (however, it may indicate a fault with the following conditions).

Not washing bedpans

1. Check that the cold-water pump is working.
2. Check that there is cold water in the tank.
3. Check sprays.
4. Check bedpan position in door cradle. (Diagram B2)

No power to control box

1. Check main power.
2. Check power to control boards.

Leaking water onto the floor

1. Check all hose fittings.
2. Check water levels and water probes.
3. Check pump seals.
4. Check door seal.

Operation, Maintenance and Installation Manual

Loading Bed Pans




Diagram B2



Side View
Diagram B3



Top View
Diagram B4

5.16 Making equipment safe after incomplete operating cycle

If an incomplete cycle has occurred, it is assumed that the load being processed has not been subjected to a complete cycle as per the set process parameters. For this reason, access to the load is restricted and will not be allowed until a full cycle has been completed.

At the end of an incomplete cycle or on powering up the device after an incomplete cycle, a fault code will display relating to the cause of the incomplete cycle (Please see section 5.13 Fault list for description of fault code). Depending on the nature of the fault an error reset may need to be completed; this can be done by following the instructions below. If the device does not allow an Error Reset to be completed, or the fault re-occurs after an Error Reset, it is recommended that the device is isolated and service is arranged to correct the fault.



Operation, Maintenance and Installation Manual

6.0 Technical Specifications

6.1 Power and Water Consumption

| POWER AND WATER CONSUMPTION DATA | | | | |
|-------------------------------------|-------------------|-------------------|-------------------|---------------------|
| Model | Avg Cycles per/Hr | Avg Cycle min/sec | Avg kWh per cycle | CW Avg Lt per cycle |
| ES915S-D / ES915P-D 20 Amp Model | 11 | 5.3 | 0.300 | 36 |
| ES935S-D / ES935P-D 10 Amp Model | 7.7 | 7.5 | 0.236 | 35.5 |




Note: Values may change due to operating and supply service conditions.

6.2 Device Specifications

| | | | | |
|------------------------------------|-------------------|---|--|-----------|
| ES-D915 Models | Electrical Rating | Volts | 240V | APPROVALS |
| | | Phase / Hz | 1 ph / 50 Hz | |
| | | Amps | 20 Amps | |
| | Heating Elements | Steam Generator | 4500W 240V 3 x 6.25A Star | |
| | | | | |
| ES-D935 Models | Electrical Rating | Volts | 240V | |
| | | Phase / Hz | 1 ph / 50 Hz | |
| | | Amps | 10 Amps | |
| | Heating Elements | Steam Generator | 2200W 240V 3 x 3.05A Star | |
| ELECTRICAL (Common components) | | | | |
| Element over temperature cut-out | | Steam Generator | Capillary bulb thermal cut out Manual reset 115°c 25A | |
| | | Water tank | Capillary bulb thermal cut out Manual reset 100°c 25A | |
| Fuse: Transformer | | Cylinder type 5 x 20mm | F3.15AL 250V | |
| Fuse F1: Power supply PCB | | Cylinder type | T315mA 250V | |
| Door Actuator | | Type LA12 Self-Locking | 24VDC | |
| Interlock Micro Switch | | SPDT | 240VAC 16A | |
| PCB HARDWARE | | | | |
| PCB (Printed Circuit Boards) | | Power supply board | CM00PS03 | |
| | | PCB module | CM00AV15 | |
| | | Interceptor PCB module | CM11IC03 | |
| | | | | |
| Environment conditions (Operating) | | Location | Indoor Use | |
| | | Altitude | <2000m | |
| | | Temperature | +10°C to +25°C | |
| | | Relative Humidity | +30% to 75% | |
| | | Mains voltage fluctuations | 240+10% | |
| | | Overvoltage | Category II | |
| | | Pollution | Degree II | |
| WASH CHAMBER | | | | |
| Loading Configurations | | 2 x Urinal Bottles and One (1) of the following: | | |
| | | Bed Pan & Lid | Commode & Lid | |
| | | Small Commode & Lid | Small Slipper Pan & Lid | |
| | | Large Slipper Pan & Lid | | |
| Chamber Specifications | | Material | 1.2mm 304 Stainless Steel | |
| | | Design Life | Minimum of 10,000 Cycles | |
| | | Usable Volume | 50 L | |
| | | Total Volume | 60 L | |




Operation, Maintenance and Installation Manual

| WASH SYSTEM | | | |
|---|---|---|--|
| Nozzles and Sprays | Flush nozzles | 1 x Rotating 180° Nozzle | |
| | Main Chamber | 2 x Fixed 60° Nozzles | |
| | Urine Bottle sprays | 2 x Fixed 60° Nozzles | |
| | Cool Down sprays | 1 x Fixed 60° Nozzle | |
| Wash pump | Horizontal multistage | 240V 1ph 3.1A | |
| Detergent | Malmest specific (5Lt) | Detergent Concentrate Detergent Pro | ARTG Class 1 |
| Detergent pump | Self-priming peristaltic pump | 240V 50Hz | |
| PLUMBING | | | |
| Soil line Connection | Type | 'S' or 'P' Trap (6mm PE) 80° | Complies to AS 2887 & AS/NZS 3500.2 |
| | Size | 100mm ID | |
| Backflow Prevention | Air Gap (RAG) | Water storage tank | AS 2845.2 |
| Cold water inlet Supply Potable Water (Final Rinse Water) | pH | 6.5-8.5pH | |
| | Chloride | <30mg/L | |
| | Water Hardness | <2.5mmol/L | |
| | Water conductivity | <850 µS/cm | |
| | Temperature | 15 – 25°C | |
| | Pressure | 100 – 400kpa | |
| | Nominal Flow Rate | 10L/Min | |
| | Solenoid Valve 1 in 2 out | 240V - GB¾ male | WRAS Certified |
| Hose – valve to Cold Water Tank | S/Steel Braided | 10mm | WRAS Certified |
| Hoses – Mains Water Supply to Inlet Valve | S/Steel Braided | G¾ Hex Nut Each End 1.5M x 10mm | Certified to AS/NZS 3499 |
| CONSTRUCTION MATERIALS | | | |
| Materials | Frame | 25 x 25 x 1.6 tube SS | |
| | Wash chamber | 304 SS grade 4 | |
| | Door | 304 SS grade 4 | |
| | External panels | 304 SS grade 4 | |
| | Water tank | 316 SS grade 2B | |
| | Steam generator | 316 SS grade 2B | |
| | Fasteners | 304/316 SS | |
| Principal Heavy Component | Main Wash Pump | 11kg | |
| TRANSPORTATION | | | |
| HANDLING & STORAGE CONDITIONS |  Fragile |  Keep away from rain |  Do not stack |
| | Temperature | -5°C to +50°C | |
| DIMENSIONS & WEIGHTS | | | |
| Weights | Nett: | Operating: | Shipping: |
| | 92 KG | 122 kg | 110 kg |
| Shipping (crated): | | 151 kg | |
| Dimensions (W x D x H) | 595 x 660 x 1300 (mm) | | 710 x 790 x 1520 (mm) |
| Floor Loading at Each Support as viewed from front of unit. * Loadings are a guide only & will vary between installations. | Front Left | Front Right | Rear Left |
| | 31 kg* | 33 kg* | 37 kg* |
| | | | Rear Right |
| | | | 21kg* |

Operation, Maintenance and Installation Manual


6.3 Electrical Diagram – 20 Amp Models (ES915 Manufactured Post 28/02/2025)





Operation, Maintenance and Installation Manual


6.4 Electrical Diagram – 10 Amp Models (ES935 Manufactured Post 28/02/2025)






Operation, Maintenance and Installation Manual

6.5 Electrical Diagram – 20 Amp Models (ES915 Manufactured Pre 28/02/2025)




Operation, Maintenance and Installation Manual

6.4 Electrical Diagram –10 Amp Models (ES935 Manufactured Pre 28/02/2025)





6.5 Hydraulic Circuit Diagram - ES915/935





Operation, Maintenance and Installation Manual

Warranty Statement

This warranty is provided, and operates in addition to, the statutory warranties Malmet (Australia) Pty Ltd ("Malmet") provides to any consumer under the Australian Consumer Law (if applicable) or by virtue of any other applicable legislation.

Subject to the following conditions, we provide, from the date of purchase, the following warranty on Malmet devices and spare parts for products manufactured by Malmet and sold in Australia:

- Functional components found within the device to be defective in workmanship or material will be repaired or replaced free of charge subject to the periods of warranty specified in the table below.
- A decision regarding whether the defective components will be repaired or replaced will be determined at the sole discretion of Malmet or its authorised agents or representatives.
- The structural warranty covers any structural components within the device, which fail to perform their intended function due to faulty manufacture or deterioration within the warranty period.
- Parts replaced in devices under warranty are warranted for the balance of the original warranty period for that device.

| Malmet Devices | |
|----------------------|-------------------------------|
| Device Components | Parts & Labour |
| Structural Guarantee | 2 Years from Date of Purchase |
| All other components | 2 Years from Date of Purchase |

| Malmet Spare Parts |
|------------------------------|
| 1 Year from Date of Purchase |

The installer is responsible for the correct installation, start up and demonstrating the operation of the product. They are also responsible for issuing the relevant certificates of compliance (these may differ from state to state).

CONDITIONS AND EXCLUSIONS

- Device must be installed and commissioned according to Malmet's instructions (outlined in Malmet Operation, Maintenance and Installation Manual) and operated to the purpose it was designed.
- Device must be serviced as instructed in the Operation, Maintenance and Installation Manuals.
- To the extent permitted by law, this warranty shall not cover damage, malfunction or failure resulting from accident, misuse or misapplication, improper or unauthorised repair, neglect or modification or use of unauthorised replacement parts or accessories, inclusive of detergent, or improper voltage. The warranty may be void if the serial number is removed or altered.
- Parts damaged in transit back to Malmet Leeton due to poor packaging could result in warranty claim being rejected in part or in full.
- Any part tampered with or which has been altered by unauthorised repairs and/or modifications will be rejected under a warranty claim to the extent permitted by law (to the extent the Australian Consumer Law applies, Malmet will assess the extent to which the tampering or unauthorised repairs contributed to the failure).
- Reasonable access must be allowed for maintenance. If any additional equipment is needed to provide access to the device, this must be provided (and paid for) by the owner.
- It is the owner's responsibility to provide safe access to the device. Malmet, or any of its authorised service agents, may refuse to perform maintenance or warranty work if access is unsafe, as determined by Malmet or any of its authorised service agents acting reasonably.
- Should a warranty claim be rejected you will be advised in writing with a full explanation of our reasons.
- Malmet have a Warranty Claim Procedure that is fair to our customers and provides an efficient system of replacement and/or repair of faulty parts. If at any time you believe we are not meeting our commitment to you please contact Malmet Head Office via email: info@malmet.com.au



Operation, Maintenance and Installation Manual

- To the extent permitted by law, no responsibility will be accepted for outside elements including, but not limited to storms, pest and vermin that may cause damage to the device.
- To the extent permitted by law, no responsibility will be accepted for damage incurred as a result of, or incidental to, electrical surges or brown outs or for any other consequential damages.
- If there is no certificate of compliance for plumbing or electrical, Malmet reserves the right to refuse service on non-compliant installations.
- To the extent permitted by law, claims for damage to contents, carpet, ceilings, foundations or any other consequential loss either direct or indirect resulting from, power spikes, incorrect operation, incorrect installation, faulty product or any other cause, are excluded.
- This warranty, and to the extent permitted by law, any warranties owed by Malmet under the Australian Consumer Law or other applicable legislation, are not transferrable and cannot be sold, assigned or transferred in any other way from the purchaser to any other person.
- To the extent permitted by law, unauthorised use of any parts that were not supplied or approved for use in the applicable device by Malmet will result in this warranty and any warranty claims applicable to that device being void.
- Warranty labour (service work) shall not include devices located outside of city metropolitan areas of Melbourne, Sydney, Adelaide, Perth and Brisbane. Costs outside these areas shall be borne by the owner. The owner shall be notified of this prior to the warranty call out.
- Warranty labour (service work) shall be performed during normal business hours (Monday – Friday 7am – 4pm), excluding public holidays.
- Warranty labour (service work) performed outside of normal business hours, shall be charged at Malmet's or its authorised representative or agent's standard after-hour labour rates.
- Warranty relating to spare parts covers parts only and does not include any associated labour costs.

To the extent permitted by law, a charge will be made for work done or a service call made where:

- There is no fault apparent with the device, as determined by Malmet or its authorised representative or agent acting reasonably.
- The defective operation of the device is due to failure of electricity or water supply.
- Defects are caused by neglect, incorrect application, abuse or by accidental damage of the device.
- An unauthorised person has attempted to repair the device.
- Harsh environmental situations including, but not limited to, water quality that may cause the water tank damage cannot be covered under this warranty



Operation, Maintenance and Installation Manual

HOW TO MAKE A CLAIM UNDER THIS WARRANTY

If you believe there is a defect in a device you have purchased from Malmet, you must notify Malmet in writing of such defect, by sending an email (**Notice of Defect**) to info@malmet.com.au prior to the expiration of the applicable warranty period set out in this warranty.

For the avoidance of doubt, Malmet must receive your Notice of Defect prior to the expiration of the warranty period.

To the extent permitted by law, Malmet will not reimburse you for any expense you incur in claiming or attempting to make a claim for repair or replacement of a component under this warranty.

Please complete details below:

| | |
|------------------------|------------------------------|
| Date Purchased: | Warranty Expiry Date: |
| | |
| Sold To: | For Service Contact: |
| | |

PROOF OF PURCHASE

Please retain your proof of purchase (receipt, invoice or commissioning certificate is accepted).

E.&O.E.

In the interest of continued product improvement, Malmet reserves the right to alter specifications without notice.

AUSTRALIAN CONSUMER LAW DISCLAIMER (APPLIES ONLY TO THE EXTENT YOU ARE A 'CONSUMER' WITHIN THE MEANING OF THE AUSTRALIAN CONSUMER LAW):

Malmet goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

Queensland and Northern NSW Distributor

EVOCARE AUSTRALIA PTY LIMITED

A.B.N. 98 078 566 604

Trading as **EVOCARE** and **L&M EQUIPMENT**

P.O. Box 145, Everton Park Qld. 4053

Ph: 07 3355 8000

Website: <http://www.evocare.com.au>

Email: sales@evocare.com.au

workshop@evocare.com.au

warehouse@evocare.com.au

accounts@evocare.com.au



Manufactured by

Malmet (Australia) Pty Ltd

ABN 95 001 717 791

www.malmet.com.au