

PELICAN POOL & SPA HOIST INSTALLATION GUIDE




Para Mobility
Making everyday living easier since 1983



INSTALLATION INSTRUCTIONS PELICAN POOL HOIST

TOOLS REQUIRED

- Hammer Drill
- 16mm masonry bit
- Hole Cleaning Pump
- Hole Cleaning Brush
- Cartridge gun – for adhesive
- Adhesive (included)
- Acetone

1. **Place hoist in position, ideally around 100 to 500mm from edge of pool** (can be further but check that lifting arm will reach into pool far enough to lower person into water). Check for obstructions when arm of hoist is rotated. Use the hoist base as a template.

2. **NB Anchors must be set into the concrete NOT into the tile bed or pavers. If there are pavers surrounding the pool a concrete pad must be laid** – contact Para Mobility for instructions.

3. **Drill 16mm diameter hole to a minimum of 100mm into concrete substrate.** (Example: 100mm **below** any tile bed).
4. **Remove dust and debris** by brushing and blowing 3 times each (If hole is wet or flooded, remove water with wet/dry vacuum). Hole may be damp but must not be flooded.
5. **Prepare all holes and socket studs** before continuing onto next step.
6. **Screw mixing nozzle onto Adhesive cartridge** (supplied) and dispense 2-3 trigger pulls of adhesive to waste until colour is grey with no streaks.
7. **Insert tip of nozzle** to bottom of hole and dispense adhesive.
8. **Fill hole** to about 2/3 full.
9. **Insert threaded stud socket**, with M12 threaded socket on top, with a rotating motion to release trapped air. Finish the top of the M12 threaded socket flush with finished floor, making sure that there is no adhesive in the threaded socket.
10. **Wait until adhesive has fully cured before loading.**
11. **Clean up** with acetone.

Earthing: As per AS3000 of the wiring rules standards – “In certain circumstances” it may be necessary to have all conductive metal objects equally potentially bonded (earthed), this includes Pool Hoists.





The pool hoist must be earthed by a qualified electrician.



**DRILL
HOLE**

**BRUSH
& BLOW**

**INJECT ABRASIVE
& INSERT STUD**



NB. Tiles can be core drilled (16mm) for a cleaner / neater finish.



IT'S ALL ABOUT ACCESS
& **FREEDOM**



Para Mobility
Making everyday living easier since 1983



Disability Discrimination Act 1992



Proud supporters of



Distributor

Evocare
www.evocare.com.au
Ph: 07 3355 8000
Email: sales@evocare.com.au