

OWNER'S MANUAL

 **Evocare**[®]
www.evocare.com.au
Ph: 07 3355 8000 Fax: 07 3355 5043
Email: sales@evocare.com.au

The BC1 Mobile Day Chair



Thank you for purchasing this Evocare BC1 Mobile Day Chair. It is strongly recommended that you read through this user manual prior to using your BC1 Mobile Day Chair for the first time. Your BC1 Mobile Day Chair should provide years of trouble-free use, provided that you adhere to these recommendations.

WE ARE HERE TO HELP YOU

To request product service and order replacement parts, please call:

Customer Service: 07 3355 8000

Spare Parts: 07 3355 8007

8.30 am to 5.00 pm Monday to Friday

Email: sales@evocare.com.au

If the product does not function properly contact Evocare Australia.

When requesting service, please provide the following information:

1. Part Code and Description
2. Date of Purchase

Address: Evocare Australia Pty Limited
260 South Pine Road, Enoggera, Qld. 4051
P.O. Box 1144, Stafford, Qld. 4053

DESCRIPTION

Product Codes:	
AA5110021 plus AA52050	In Ink Vinyl Upholstery
AA5110022 plus AA52050	In Latte Vinyl Upholstery
AA5110040 plus AA52050	In Crypton Birch Upholstery
AA5110041 plus AA52050	In Crypton Jasper Upholstery
AA5110042 plus AA52050	In Crypton Pimento Upholstery
AA52050	Kit to make mobile. Includes Castors, Push Handle and Slide out Foot Rest
AA50904	Tray
AA50950	Tray attaching Kit

Features

- With 100mm castors, 2 x swivel, 2 x total lock (lockable castors at rear as standard).
- A one piece ergonomically designed shell.
- Push handle
- Slide out footrest
- The frame is fully welded
- 10 year structural warranty
- ARTG number: **218697**

Optional

- Tray and tray storage
- 125mm castors

The BC1 Mobile Chair Specifications	
Overall Chair Width:	730mm
Overall Depth:	740mm
Overall Height:	990mm to 1140mm
Seat Width:	520mm
Seat Depth:	470mm
Seat Height:	360mm to 510mm
Back Height:	720mm
Safe Working Load (SWL):	150kg

SAFETY

SAFETY:

These instructions are guidelines only, however it is recommended that you read and understand these instructions before using the BC2 Mobile Day Chair. Failure to do so may cause injury or damage.

Maximum user weight 180kg



WARNING All the BC2 Mobile Day Chair wheels must be firmly in contact with the floor and brakes must be applied when sitting and standing.



WARNING The moving actions of chairs can be dangerous to young children. Always ensure children are supervised when in the vicinity of the product in use



WARNING Do not place the chair over power leads



WARNING Do not permit more than one person on the BC2 Mobile Day Chair at one time



WARNING Check regularly for wear and tear.

CLEANING & MAINTENANCE

CLEANING

- Vinyl covers should be wiped over using a mild solution of soap and warm water.
- The life of the fabric and upholstery can be increase by regular vacuuming.
- All leather chairs should be treated regularly with a leather conditioner.
- Standard fabrics can be treated with any protector to reduce the risk or permanent staining.

MAINTENANCE

- Check castors twice per year to ensure free movement.
- Check sliding footrest twice per year to ensure free movement.
- Lubricate footrest slider once per year to ensure free movement.

WARRANTY

1. Evocare provides consumers with the following warranty in relation to this Product, in addition to complying with the requirements of any relevant legislation, including the *Competition and Consumer Act 2010* (Cth) in Australia and the *Consumer Guarantees Act 1993* in New Zealand (the Acts), except where a New Zealand consumer acquires the Product for the purposes of a business.
2. In this warranty, we have used the following definitions:
 - (a) **Form** means a warranty claim form issued by Evocare in respect of Products;
 - (b) **Evocare or our** means Evocare Australia Pty Limited, ABN 98 078 566 604 of Level 19, 144 Edward Street, Brisbane 4000;
 - (c) **Products** means the goods manufactured for Evocare (including products manufactured by its contract manufacturers);
 - (d) **Material** means a material or component used by Evocare in the manufacture of the Products;
 - (e) **Retailer** means the authorised dealer of Products from whom the Product was purchased;
 - (f) **Warranty Period** means the warranty period commencing from the date of purchase of the relevant Product. If the Product or part (as the case may be) is repaired or replaced, there will be no extension to the original warranty period; and
 - (g) **Workmanship** means the handling, assembly and manufacturing processes performed by or on behalf of Evocare in order to manufacture the Products.
3. Evocare warrants that for the duration of the Warranty Period, all Products will be free of faults arising from defects in Workmanship or Materials, on the terms and conditions set out in this warranty.
4. Evocare undertakes that if during the Warranty Period any Product, or any part of a Product, has failed to operate correctly due to faulty Workmanship or defective Material, it will repair or replace the Product or part (as the case may be) free of charge provided that the following procedures are met :
 - (a) The consumer must contact the Retailer or Evocare upon becoming aware of any defect to a Product upon which the consumer will be required to complete a Form and provide satisfactory proof of purchase.
 - (b) Evocare will review the Form to determine whether there is a defect, and if so Evocare agrees to (at its option) repair, replace or supply equivalent goods, or pay the cost of any of those remedies to the consumer.
 - (c) If Evocare requests the return of the applicable Product or part, the consumer will be responsible for the collection and freight costs of returning that Product or part to Evocare.
 - (d) The consumer acknowledges that it is also responsible for the freight costs to deliver any new Product or replacement part to it.
5. This warranty is in addition to any non-excludable legal rights or remedies conferred on the consumer under any applicable Act and any similar laws. To the extent permitted by law, Evocare's liability for any non-excludable condition or warranty is limited to rectifying any defect at its option, as set out in paragraph 4(b).
6. Subject to the requirements of any applicable Act or legislation and to the extent permitted by law, no liability (whether expressed or implied) of any nature whatsoever, is accepted by Evocare for any consequential loss, damage or injury arising as a result of any fault in the Products.
7. This warranty does not extend to damage to Products which occurs during transit or transportation, or which is caused by any abuse, accident or improper installation, connection, use, adjustment or repair or use of goods otherwise than in accordance with instructions issued by Evocare.
8. The warranty on Products is waived if any addition or attachment to the Products do not have Evocare's approval or are not sold as Evocare products. The Products are designed to perform specific tasks under established test loads and unauthorised attachments may produce stresses for which the design is not appropriate.
9. The following applies to consumers who purchased a relevant Product in Australia:

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

WARRANTY, continued

The manufacturer warranty for BC1 Mobile Day Chair is 10 years structural.

NOTE: This does not include staining or accidental damage from misuse or neglect.

Evocare must inspect the product within the warranty period for any claims lodged.

Evocare warranty does not cover product that has been modified from its original form.

For more information about your BC1 Mobile Day Chair, warranty or to order parts call 07 3355 8007.

E- mail: sales@evocare.com.au
Website www.evocare.com.au