



**MALMET (AUSTRALIA)
PTY LTD**

Head Office and Customer Service
ABN 95 001 717 791

9-11 McKay Avenue
PO Box 373
Leeton NSW 2705

Phone: +61 2 6953 7677
Email: info@malmet.com.au

Bedpan/Urinal Bottle and Utensil/Bowl Washer Disinfector

Models WDS1; WDS3



Operation, Maintenance and Installation Manual

**Note: Due to Malmet's Policy of continuous product improvement;
design and technical specifications are subject to change without notice.**

Serial Number:	Supplied to:
Date Installed:	Installed by:

It is important that the name from whom you purchased your device and the name of the installer are recorded above. The installer is responsible for the correct installation, start up and demonstrating the operation of this device. They are also responsible for issuing relevant certificates of compliance (these may differ from state to state).



Table of Contents

Foreword	2
Certifications and Compliances.....	2
Quality Policy.....	2
Important Warranty Reminder	2
Malmet Head Office and Factory Contact Details	2
Safety Instructions – Warnings	3
Intended Use	4
1.0 Design Parameters.....	4
1.1 Operating Cycles	5
1.2 Chemical Dosing System.....	6
1.3 Detergent.....	6
1.4 Device Features.....	8
1.5 Control LCD Display Features	8
1.6 Operating Features	9
2.0 Installation and Commissioning.....	10
2.1 Installation	10
2.2 Service Connections.....	11
2.3 Plumbing Connections	11
2.4 Venting.....	13
2.5 Electrical.....	13
2.6 Service Connection Points	14
2.7 Device Dimensions.....	15
2.8 Commissioning (To only be completed by qualified persons)	16
3.0 Loading and Operation.....	17
3.1 Urinal Bottle / Bedpan – Loading (Refer to laminated quick reference instructions)	17
3.2 Urinal Bottle / Bedpan – Operating Cycle (Refer to laminated quick reference instructions)	18
3.3 Urinal Bottle Only – Loading (Refer to laminated quick reference instructions).....	19
3.4 Urinal Bottle Only – Operating Cycle (Refer to laminated quick reference instructions).....	19
3.5 Bowl / Utensil – Loading (Refer to laminated quick reference instructions).....	20
3.6 Bowl / Utensil – Operating Cycle (Refer to laminated quick reference instructions).....	21
4.0 Cycle of Operation	22
4.1 Release of Processed items	27
5.0 Maintenance	27
5.1 Daily Maintenance (Operator or Maintenance Technician)	27
5.2 Bi-Monthly Maintenance (Maintenance Technician).....	27
5.3 Recommended Preventative Maintenance Schedule (to be performed by qualified maintenance personnel)	28
5.4 Post Maintenance or Repair Safety Checks.....	29
5.5 Door Safety Test	29
5.6 Cycle Operation Check.....	29
5.7 Post Maintenance and Repair of Electrical Systems	29
5.8 Replacement of Safety Devices	29
5.9 Validation.....	30
5.10 Requalification.....	30
5.11 Device faults.....	31
5.12 Fault Indication	31
5.13 Fault Codes.....	31
5.15 Making equipment safe after incomplete operating cycle	34
5.16 Purge Reset (Operator).....	34
5.17 Purge Detergent Line (Operator)	34
5.18 Manual Purge Reset (Maintenance).....	35
5.19 Manual Purge Detergent Line (Maintenance).....	35
5.20 Safe Access to the Load or Chamber in Fault Condition.....	36
6.0 Technical Data	37
6.1 Power and Water Consumption	37
6.2 Device Specifications	38
6.3 Wiring Diagram (1ph devices) Manufactured after June 2020	40
6.3 Wiring Diagram (1ph devices) Manufactured before June 2020	41
6.4 Wiring Diagram (3ph devices) Manufactured after June 2020	42
6.4 Wiring Diagram (3ph devices) Manufactured before June 2020	43
Warranty Statement	44



Foreword

In order to obtain maximum life and efficiency from your Malmet Bedpan/Urinal Bottle and Utensil/Bowl Washer Disinfectant and to aid in the safe operation of the device, please read and understand this manual thoroughly and follow all instructions before operating.

This device is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge.

Children should be supervised to ensure that they do not play with the device.

The specifications supplied in this manual were in effect at time of publication, however due to Malmet (Australia)'s policy of continuous improvement, changes to these specifications may be made at any time without notice on the part of Malmet (Australia) Pty Ltd.

Certifications and Compliances

ARTG Identifier: 232450

Electrical Safety Cert: Cert No. CS10462N IEC61010-2-040 in conjunction with IEC61010-1

EMC Compliance: Cert No: M1910013

Watermark Cert: Cert No: WMKA21156

Quality Policy

Malmet's quality management system is certified to ISO 13485 and ISO 9001 and guarantees the quality of this device.

Important Warranty Reminder

Should you have any problems with your device, contact the company from whom you purchased it, or Malmet (Australia) Pty Ltd.

It is important that the name from whom you purchased your device and the name of the installer are recorded on the front page of this manual. The installer is responsible for the correct installation, start up and demonstrating the operation of this device. They are also responsible for issuing relevant certificates of compliance (these may differ from state to state).

Malmet Head Office and Factory Contact Details

Malmet (Australia) Pty Ltd

9-11 McKay Avenue
PO Box 373
LEETON NSW 2705

Telephone: +61 2 6953 7677

E-mail: info@malmet.com.au

Website: www.malmet.com.au

Safety Instructions – Warnings

Please read and understand this manual before using this device, if this device is used in a manner not specified by the manufacturer protection by the device may be impaired.

Please refer to this manual for information wherever this warning symbol is displayed –



Be aware of 240V / 415V Voltage



Disconnect power when servicing



Mains power ISO switch must be in an accessible position easily reached and not obstructed so device can be isolated from mains power during service



For the safe use of this device the responsible body should ensure that all operators are adequately trained to operate and maintain the device in its safe use.



Only process items listed as per the design parameters in this manual



Be aware of steam discharge



Utensils and racks are hot to handle



Safety gloves and goggles must be worn when changing detergent



Safety clothing with reflective tape can activate the hands-free sensor when device is in standby mode



Be aware of hot surfaces, pipes and hoses from steam and hot water



Install temperature probes and element over temperature thermal cut-outs correctly



Plumbing service connection must comply with AS/NZS 3500.1 and AS/NZS 3500.2



Not suitable for use in the presence of flammable anaesthetic mixtures with air or nitrous oxide and mode of operation as continuous



Intended Use

The Malmet Bedpan / Urinal Bottle and utensil/Bowl Washer Disinfector is intended to be used by health service organisation personnel for the cleaning and disinfection of medical devices during an automatic cycle of the types intended to be re-used as listed in the design parameters section of this manual.

1.0 Design Parameters

The Malmet Bedpan/Urinal Bottle and Utensil/Bowl Washer Disinfector has been designed within the following parameters:

- a) To empty, clean and disinfect devices of the type listed below intended to be re-used and that have been soiled by human excreta. *
- b) To clean and disinfect during each automatic cycle a load comprised of the following re-useable medical devices.
 - i) Standard size bedpans
 - ii) Standard male and female urinal bottles
 - iii) 300mm/255mm/220mm/160mm Kidney dishes
 - iv) 345mm/305mm/240mm/210mm/185mm/140mm/110mm/100mm/80mm bowls, dishes and receivers.
- c) Two bedpans with lids and four urinal bottles can be emptied, cleaned and disinfected during an automatic cycle.
- d) The cradle is designed to ensure that re-useable medical devices are not dislodged during the cleaning cycle. The contents are emptied during door closure.
- e) The chamber and door are self-cleaning and do not permit water or soil to remain after a properly completed cycle. Steam disinfecting ensures all internal surfaces are totally clean and safe.
- f) The bedpan flush cycle:
 - i) Removes the heavy soil
 - ii) Clears the trap
- g) Disinfection cycle is factory pre-set to 90° for 1 minute, to achieve an A₀ of 600 in accordance with AS 5369:2023, Section 6.3, Table 6.1.

*The efficacy of disinfection can be impaired if soil removal is incomplete before the start of the disinfection process. Users should be aware that some medical devices might require pre-treatment such as soaking or brushing before processing. Please follow the manufacturer's instructions when processing re-useable medical devices.

**High disinfection temperatures may affect some heat sensitive reusable medical devices; ensure items are disinfected as per the manufacturer's instructions.

*Non-flushable items such as Feminine Hygiene products, Incontinence or Absorbent pads & diapers, Condoms, Dental Floss, Paper Towel, Cigarette butts, medications or chemicals cannot be disposed of in Malmet Washer Disinfectors and should be disposed as per facility &/or local waste management requirements.



1.1 Operating Cycles

Three available operating (cleaning and disinfection)



Urinal Bottle / Bedpan

Load capacity: 2 x small slipper pans, 2 x large slipper pans, 2 x standard bedpans and 4 x standard male and female urinal bottles

Cycle:

- | | |
|----------------------------------|---------------|
| 1. Flush | 8 to 12 sec |
| 2. Cold water wash | 25 to 35 sec |
| 3. Hot water wash with detergent | 10 to 15 sec |
| 4. Hot water rinse | 25 to 35 sec |
| 5. Disinfection | 1 min at 90°C |
| 6. Cool down | 20 sec |



Urinal Bottle only

Load capacity: 4 x Standard male and female urinal bottles

Cycle:

- | | |
|----------------------------------|---------------|
| 1. Cold water wash | 40 to 50 sec |
| 2. Hot water wash with detergent | 40 to 50 sec |
| 3. Hot water rinse | 25 to 35 sec |
| 4. Disinfection | 1 min at 90°C |
| 5. Cool down | 20 sec |



Bowl / Utensil

Load capacity: 2 x large bowls (345Dia to 305Dia), 3 x medium bowls (240Dia to 210Dia), 6 kidney dishes (sizes 300 x 50, 255 x 50 and 220 x 43 – held in kidney rack holder)

Smaller kidney dishes fit into a basket with a smaller open-ended trays and utensils.

Cycle:

- | | |
|----------------------------------|---------------|
| 1. Cold water wash | 40 to 50 sec |
| 2. Hot water wash with detergent | 40 to 50 sec |
| 3. Hot water rinse | 25 to 35 sec |
| 4. Disinfection | 1 min at 90°C |
| 5. Cool down | 20 sec |



Note: These times do not include filling and heating.

In accordance and complies with:

Cleaning Efficacy
AS 5369, ISO 15883-1, ISO 15883-5

Thermal Disinfection
AS 5369, ISO 15883-1



1.2 Chemical Dosing System

The chemical dosing system ensures the correct amount of process chemicals are delivered at the correct stage during each operating cycle. Chemical addition is controlled by the automatic cycle controller and will display a fault if the incorrect amount of chemical is injected. The dosing system is designed for use with Malmet Chemical; other chemicals may impair the operation of the system.

1.3 Detergent

Malmet Washer Disinfectors are designed and tested using Malmet specific detergent. To ensure intended performance and trouble-free operation use only Malmet branded detergent, failure to use correct detergent may void warranty.

The detergent system delivers 30ml of concentrated detergent during the Hot Water Wash stage.

Changing the detergent bottle



WARNING!



SAFETY GLOVES AND GOGGLES MUST BE WORN WHEN CHANGING DETERGENT AND DISPOSING OF EMPTY CONTAINERS

Current Safety Data Sheet for Malmet detergent is available in .PDF format from Malmet's website. www.malmet.com.au

The 5-litre detergent container is accessed by opening the bottom door. Only use Malmet approved detergent (See technical data for detergent details)

1. Pull handle on detergent door and open
2. Unscrew cap and pull out with suction hose (let hose hang in detergent chamber)
3. Remove empty bottle and replace with full bottle

Note: Leave cap on new bottle until in position

4. Remove cap on new bottle and fit existing hose and cap, make sure suction hose is at bottom of bottle
5. Close detergent door
6. Restart device operation as normal




Disposal of Empty Detergent Bottles

1. When handling empty containers treat as though container is full and wear appropriate PPE as per MSD requirements
2. Triple rinse container and puncture base to render unusable
3. Dispose of containers as per facility chemical waste handling procedures, and/or local government requirements



Detergent Warnings

If insufficient detergent is available to deliver the required quantity, the device will display a fault 302. Access to the load will be restricted until the condition is rectified.

To rectify this fault, install a new detergent bottle as per the procedure in this manual “Changing the Detergent Bottle”.



Once a full bottle has been installed the detergent line will require purging

1. Turn power OFF  at control LCD display
2. Turn power ON  at control LCD display


On powering back on the display will show TEST 302



3. The detergent line will automatically purge if necessary for a predetermined time or until the detergent is detected
4. Once detergent is detected the fault will clear and the previous cycle will be restarted

If the detergent fault is not cleared check detergent level and repeat the purge process again, if the detergent fault continues and is unable to be cleared contact maintenance or service provider

1.4 Device Features











1.5 Control LCD Display Features







1.6 Operating Features

	POWER	On/Off Standby
	HANDS FREE SENSOR	For hands free opening, closing door and starting of cycle. To operate, hold hand in front of sensor
	URINAL BOTTLE ONLY	Select and start for urinal bottle only cycle
	BOWL / UTENSILS	Select and start for bowls / kidney dishes etc
	URINAL BOTTLE / BEDPAN	Select and start for urinal bottles / bedpans
	DOOR	To manually open and close door, if pressed whilst door is closing at the start of a cycle the cycle will be cancelled and the door opened
	SCROLL	Menu scroll button SERVICE ONLY
	SELECT	Menu select and enter button SERVICE ONLY Purge Reset: When “purge” is flashing on display, if pressed this will start the purge cycle after a fault.

Hands free operation

	The hands-free option can open and close the door without touching the device
	Adjusting the rack position and loading the device accordingly (refer to 3.0 Loading and Operation)
	The hands-free option can be used to start the operating cycle

Manual operation of the door

	Press the door button to open the door
	Adjust the rack position and load the device accordingly (refer to 3.0 Loading and Operation)
	Press the door button to start the operating cycle

2.0 Installation and Commissioning

2.1 Installation

To avoid problems with this device, these Installation Guidelines should be followed.


Installations must be carried out by a qualified and licenced tradesperson.

The device must be affixed to the building structure to ensure mechanical stability prior to use. The device can be affixed to the building structure via the holes provided in the device feet, or use of an appropriately rated appliance safety strap.

Prior to installation of the device, services as noted are to be provided by the facility. It is not the responsibility of Malmet to provide these service connections.

Risk Assessment

It is recommended a risk assessment is conducted by the user both prior to and after installation and any risks identified mitigated to an acceptable level using the hierarchy of control;



<https://commons.wikimedia.org/w/index.php?curid=55610678>

Handling



Weights of Device: Net: 140 kg Shipping: 165 kg Shipping with crating: 217 kg.

- Handling of the device to installation site must be with fork lift or hand pallet truck.
- Before unpacking device inspect carton for any damage relating to forklift forks and damage relating to device falling over or for evidence of top loading
- After unpacking the device, inspect all external panels for damage.
- Remove the 4 screws holding the device to the pallet.
- Follow your internal manual handling guidelines to manoeuvre the device off the pallet. The device can then be placed into position by fork lift or hand pallet truck.



Disposal of Packaging

- Please dispose of packaging as per facility procedures or local government requirements.



Disposal of Medical Device

- Please dispose of medical device as per state environmental regulatory requirements.


2.1.1 Positioning the WDS

Model	Placement	Access Required	Device Dimensions		
			Height	Width	Depth (mm)
WDS1; WDS3	Freestanding	Both Sides	1695	610	625

Please allow sufficient room for servicing purposes. Minimum recommended clearance of 150mm on either side and 150mm at the rear of the device. Clearances below the recommended measurement can be achieved however additional charges may be incurred for servicing and maintenance of the equipment.

New Buildings

- Service connections are usually pre-placed after planning and consultation with all interested parties. Installation is by connection to the services provided.
- As the soil line (sewerage outlet) is the least flexible of all the connections, this usually influences the decision as to where to place the WDS. If an existing soil line can be utilised this will represent a cost saving.
- The WDS is supplied as either a 'S' or 'P' Trap as nominated by the Purchaser. The 'S' Trap connects through the floor and the 'P' Trap connects through the back wall. The trap section is easily removed if the wrong trap has been ordered. Refer to Diagram C2 for trap connections.

 Potential electromagnetic or other interference between other EQUIPMENT and devices can possibly affect the infra-red hands-free operation sensor. It is advisable to check all the equipment and devices in the intended installation area that have infra-red operation. **Electromagnetic interference can be prevented by installing the device in non-patient areas of hospital (or similar).**

2.2 Service Connections

MODEL	HOT WATER	COLD WATER	SOIL LINE	ELECTRICAL
WDS1 (1ph)	Solenoid valve GB ¾ Male	Solenoid valve GB ¾ Male	100mm 'S' or 'P' Trap Hot water discharge 50-55°C	240V 1 Phase 20 Amps 50 Hertz
WDS3 (3ph)	Solenoid valve GB ¾ Male	Solenoid valve GB ¾ Male	100mm 'S' or 'P' Trap Hot water discharge 50-55°C	415V 3 Phase 20 Amps 50 Hertz

2.3 Plumbing Connections



These installation guidelines must be followed to ensure the device will operate as intended.

Installations must be carried out by a qualified, licensed tradesperson.



Water Supply – Note: Water Supply connections must comply with AS/NZS3500.1 and be Watermark Certified.

HOT AND COLD-WATER CONNECTIONS ARE REQUIRED. The device can be connected to any potable mains pressure hot and cold-water supply as a back-flow prevention air gap is incorporated in the design of the water tank. Complies to AS 2845.2

Cold Water – Note: See device specifications table for water quality requirements

Pressure: 100-350kPa

Temperature: 15-25°C

Nominal Flow: 10L/Min

Connection: GB ¾" BSP

Hot Water – Note: See device specifications table for water quality requirements

Pressure: 100-350kPa

Temperature: 55-65°C


Nominal Flow: 10L/Min

Connection: GB ¾" BSP

- The water supply is to be connected to an isolating valve or cistern stopcock placed approximately 1200mm from the floor to the right-hand side of the device (preferably not behind the device).
- If inlet water flow pressure is higher than 350 Kpa an inline pressure reducing valve should be fitted

Note: Old hose sets should not be re-used; it is recommended new compliant hose sets are used on installation.

Waste / Soil Line Connection



For reasonable connection working space allow 150mm from rear of device to wall.

The device is 610mm wide and the centre of the trap is 305mm from each side.

FREE STANDING

'S' TRAP PIPE POSITIONING

For optimal installation, the center of the soil line designated for the 'S' Trap should be positioned at an approximate distance of 227.5mm from the back wall.

Additionally, for adequate side service access, a spacing of 455mm from the soil line to the nearest side wall is advisable. In cases where spatial constraints limit adherence to these recommendations, priority should be given to maximizing space on the right-hand side (as viewed from the front of the device), maintaining a minimum clearance of 150mm. This arrangement is recommended to facilitate easier servicing of the steam tank element and probe, minimizing potential difficulties.

FREE STANDING

'P' TRAP PIPE POSITIONING

The centre of the soil line to receive the 'P' Trap should be approximately 410mm from the floor when the device is positioned 150mm from the wall. Because this pipe is graded to 5° this measurement will vary as the device is installed closer or further away from the back wall.

- Hot water Discharge temp 50° – 55°C. Soil line must comply to AS/NZS 3500.2
- Soil line connection is by a pan collar or other preferred method. If the belled end on the polyethylene moulded trap is not required it can easily be cut off to provide a straight pipe connection.
- Level the device by using the flanged screw in legs and, if possible, maintain approx. 100mm floor clearance for ease of floor cleaning. Malmet recommends affixing some of the leg flanges to the floor via stainless self-tapping screws to prevent sideways movements and damage to services and soil line connections.
- The soil line should protrude from the floor or wall at a minimum of 100mm.

2.4 Venting

No external vent pipe work is required as the device is designed to condensate all visible steam internally

2.5 Electrical



These installation guidelines must be followed to ensure the device will operate as intended and must be carried out by a qualified, licensed electrical tradesperson. The device must be installed and serviced to national wiring rules AS/NZS 3000.

Model: WDS 1ph 240V 50Hz 1 Ph 20Amp


WDS 3ph 415V 50Hz 3 Ph 20Amp




Overcurrent Protection device: A 20A circuit breaker or fuse must be installed in the building installation for the device.

- Devices are supplied with 1700mm power supply cord extending from the rear top right-hand side of the device for hard wiring into Mains ISO switch. Position switch approx. 1500mm above floor level.
- An ISO switch or circuit breaker must be included in the installation (not supplied by Malmet).
- Mains power ISO switch must be suitably located and easily reached, approximately 1500mm above floor level adjacent to device. The 1.7m mains power lead exits the device approximately 1600mm above floor level on the right-hand side of the device.
- Mains ISO switch must be marked as the disconnecting device for the equipment.
- Mains ISO switch must have contact disconnection of all poles to provide full disconnection.
- Mains ISO switch must be suitably placed so device will not obstruct safe access to the isolator once installed.
- Ensure isolator is suitably placed so device will not obstruct safe access to the isolator once installed.
- Device must be earthed.
- Malmet recommends having a 30mA RCD in the mains supply fixed wiring.
- If the supply cord is damaged, it shall be replaced by the manufacturer or its service agent or similarly qualified person in order to avoid any potential hazard.

2.6 Service Connection Points




Dimension Tolerance $\pm 5\text{mm}$



DETAIL A
SCALE 1 : 5




2.7 Device Dimensions



Dimension Tolerance $\pm 5\text{mm}$

FRONT VIEW

RIGHT VIEW



BACK VIEW



2.8 Commissioning (To only be completed by qualified persons)

- a) Before switching on the device make sure the DEVICE IS LEVEL and WATER TAPS ARE ON. Check that the DRAIN WASTE is connected.
- b) Turn on the power at the isolation switch and press the standby button on the front display to power on. The LCD display will illuminate and go to standby mode.
- c) Check that the water tank is filling with cold water and that it fills to Level 3. Ensure the lid is put back on water tank.
- d) Check that the steam generator tank has filled to the full Level 3. The steam generator will preheat to 85°C in standby mode or during cycle.
- e) Check that the 5-litre detergent bottle has been fitted. From the menu, run the DETERG/START to purge detergent through the line to flow sensor. This may need to be done twice if the line is empty on start-up (see 5.19).
- f) Flush approximately 1 litre of water down the steam generator tank overflow pipe. This will fill the 'S' Trap at the hose junction and prevent steam coming back up into the water tank.
- g) From the Menu, run Purge Reset three times (see 5.18). Check for any water leaks.
- h) Open the door and turn power off at the control. Check racks for freedom of movement. Check that flush nozzles rotate freely. Check that all spray nozzles are tight. Turn power back on and close door by pressing the door button.
- i) Select the BEDPAN/URINE cycle and let it run through the cycle. Repeat the same procedure for the URINE and BOWL cycles.
- j) Replace all covers.

Note: DO NOT USE THE DEVICE WITHOUT THE WATER SUPPLY TURNED ON

3.0 Loading and Operation



For the safe use of this device the responsible body should ensure that all operators are adequately trained to operate and maintain the device in its safe use.

3.1 Urinal Bottle / Bedpan – Loading (Refer to laminated quick reference instructions)

1. To open door, use  press 

2. Place racks into the correct positions

Lift bedpan rack up and lock into position



Check that bowl rack is in the up position

3. Place bedpan lids into back of bowl rack



4. Place first bedpan or large/small slipper pan onto top rack










5. Place second bedpan or large/ small slipper pan onto bottom rack



6. Place urinal bottles into rack



3.2 Urinal Bottle / Bedpan – Operating Cycle (Refer to laminated quick reference instructions)

1. Open door, using  or press 
2. Set rack position and load items to be processed as shown in 3.1
3. Select URINE / BEDPAN cycle  to close door & start cycle or
4. Close door, using either the  or press 
5. When the door closes, to start cycle use the  or press  within 10 Seconds whilst “START” is flashing



Cycles stages will be shown on LCD display

1. Flush
 2. Cold water wash
 3. Hot water wash with detergent
 4. Hot water rinse
 5. Disinfection
 6. Cool down
 7. Completed
8. When the cycle is completed the display will flash *COMPLETED OPEN DOOR*


Open door, use the  or press 

NOTE: If temperature inside the chamber is too high the display will flash, *WAITING COOLING*

When the temperature has dropped to a safe level the display will flash, *OPEN DOOR*


Open door, use  or press 

To run the same cycle, repeat steps 1 to 3

If you don't need to run another cycle, close door by pressing 



NOTE:



- In an emergency, pressing the Power button  on the top left-hand corner of the display will power off the device.
- On powering up after a cancelled cycle the device will display a Fault code 904 – Power Interrupted (see note below)



NOTE:

- On powering up device will prompt to complete a “Purge Reset”, if the condition that caused the fault has been rectified press “Select” to start the “Purge Reset”. On completion of the “Purge Reset” device will restart previous cycle and will not allow access to the load until cycle has been successfully completed.

3.3 Urinal Bottle Only – Loading (Refer to laminated quick reference instructions)

1. Open door, use  or press 
2. Place racks into correct positions

Bedpan rack must be in the down position










Check that the bowl rack is in the up position

3. Place urinal bottles into rack





3.4 Urinal Bottle Only – Operating Cycle (Refer to laminated quick reference instructions)

1. Open door, using  or press 
2. Set rack position and load items to be processed as shown in 3.3
3. Select URINE ONLY cycle  to close door & start cycle or
4. Close door, using either  or press 
5. When the door closes, to start the cycle use the  or press  within 10 Seconds whilst "START" is flashing

Cycles will be shown on LCD display


1. Cold water wash
 2. Hot water wash with detergent
 3. Hot water rinse
 4. Disinfection
 5. Cool down
 6. Completed
7. When the cycle is completed the display will flash *COMPLETED OPEN DOOR*

Open door, use  or press 

NOTE: If temperature inside the chamber is too high the display will flash *WAITING COOLING*

When the temperature has dropped to a safe level the display will flash *OPEN DOOR*

To open door, use



or press




To run the same cycle, repeat steps 1 to 3

If you don't need to run another cycle close door by pressing



NOTE:

- In an emergency, pressing the Power button  on the top left-hand corner of the display will power off the device.
- On powering up after a cancelled cycle the device will display a Fault code 904 – Power Interrupted (see note below)




NOTE:

- On powering up device will prompt to complete a “Purge Reset”, if the condition that caused the fault has been rectified press “Select” to start the “Purge Reset”. On completion of the “Purge Reset” device will restart previous cycle and will not allow access to the load until cycle has been successfully completed.

3.5 Bowl / Utensil – Loading (Refer to laminated quick reference instructions)

1. Open door, use



or press



2. Place racks into the correct positions

Place bedpan rack
into the down
position



Place the bowl
rack in the down
position

Examples of bowls and utensils placed onto racks

Large and medium size bowls



Medium size kidney dishes
(Optional accessory Part # 92-5033)

Kidney dish rack placed onto bowl rack



Load kidney dishes into rack










Smaller items can be placed into a basket
(Optional accessory Part # 92-5235)

Basket placed onto bowl rack





3.6 Bowl / Utensil – Operating Cycle (Refer to laminated quick reference instructions)

1. Open door, using  or press 
2. Set rack position and load items to be processed as shown in 3.5
3. Select BOWL / UTENSIL cycle  to close door & start cycle or
4. Close door, using either the  or press 
5. When the door closes, to start the cycle use the  or press  within 10 Seconds whilst "START" is flashing

Cold water wash



1. Hot water wash with detergent
2. Hot water rinse
3. Disinfection
4. Cool down
5. Complete

When the cycle is completed the display will flash *COMPLETED OPEN DOOR*


5. Open door, use  or press 

NOTE: If temperature inside the chamber is too high the display will flash *WAITING COOLING*

When the temperature has dropped to a safe level the display will flash *OPEN DOOR*


6. To Open door, use  or press 

7. To run the same cycle, repeat steps 1 to 3

8. If you don't need to run another cycle close door by pressing 



NOTE:

In an emergency, pressing the Power button  on the top left-hand corner of the display will power off the device

- On powering up after a cancelled cycle the device will display a Fault code 904 – Power Interrupted (see note below)




NOTE:

- On powering up device will prompt to complete a "Purge Reset", if the condition that caused the fault has been rectified press "Select" to start the "Purge Reset". On completion of the "Purge Reset" device will restart previous cycle and will not allow access to the load until cycle has been successfully completed

4.0 Cycle of Operation

1. Press POWER ON button; display shows:

POWER




After short delay display shows alternating

**OPEN
DOOR**




2. Open door; display shows:

OPENING




3. When door is open; display shows:

**LOAD
ITEMS**



After a short delay, display shows the cycle selection based on the rack positions (Urine, Bedpan & Urine or Bowl)

**SELECT
CYCLE**



If the wrong cycle has been selected for the position of the racks; an audible alarm will sound and display briefly shows: Check rack position and cycle selected match load type being processed.

**WRONG
RACK**

Then


**SELECT
CYCLE**






4. After select cycle; display shows:

CLOSING



5. Once door is closed using Hands Free or door button; display flashes 'START' for 10 seconds.

START




If cycle isn't initiated within this period operation returns to Step 1

6. When door closes via cycle select or cycle initiated in 10 second period; display shows:

BEDPAN/URINE BOTTLE
CYCLE 1: COLD FLUSH

URINE ONLY OR BOWL


*COLD
FILL 1*



NOT USED

When filled to correct level; display shows:

FLUSH 1




NOT USED

Display shows:


CYCLE 2: COLD WASH

CYCLE 1: COLD WASH

*COLD
FILL 2*




*COLD
FILL 1*




When filled to correct level; display shows:

*COLD
WASH 2*



*COLD
WASH 1*






BEDPAN/URINE BOTTLE
CYCLE 3: HOT WASH PLUS DETERGENT


Display shows:

**HOT+DET
FILL 3**




**HOT+DET
FILL 2**

URINE ONLY OR BOWL
CYCLE 2: HOT WASH PLUS DETERGENT




When filled to correct level; display shows:

**HEATING
WATER 3**




**HEATING
WATER 2**




When water reaches temp; display shows:

**HOT+DET
WASH 3**




**HOT+DET
WASH 2**



CYCLE 4: HOT RINSE


Display shows:

**HOT
FILL 4**




**HOT
FILL 3**

CYCLE 3: HOT RINSE




When filled to correct level; display shows:

**HEATING
WATER 4**



**HEATING
WATER 3**




When water reaches temp; display shows:

**HOT
RINSE 4**



**HOT
RINSE 3**





7. Disinfection; display shows:

TEMP 50° to 93°C
Steam heating chamber to 93°C
When chamber reaches 93°C
Disinfection for 1 min @ 93°C

BEDPAN/URINE BOTTLE

ON
BEDPAN & URINE BOTTLE
TEMP 50°C **88**
OPERATING MINUTES TO END

URINE ONLY OR BOWL

ON
URINE BOTTLE BOWL
TEMP 50°C **88**
OPERATING MINUTES TO END

TEMP 93°C

ON
BEDPAN & URINE BOTTLE
DISINFECTION
TEMP 93°C **88**
OPERATING MINUTES TO END

ON
URINE BOTTLE BOWL
DISINFECTION
TEMP 93°C **88**
OPERATING MINUTES TO END

When disinfection completes; display briefly shows:

SUCCESS

ON
BEDPAN & URINE BOTTLE
DISINFECTION
SUCCESS **88**
OPERATING MINUTES TO END

ON
URINE BOTTLE BOWL
DISINFECTION
SUCCESS **88**
OPERATING MINUTES TO END

8. Cool Rinse; display shows:

**COOL
RINSE**

ON
BEDPAN & URINE BOTTLE
COOL-RINSE **88**
OPERATING MINUTES TO END

ON
URINE BOTTLE BOWL
COOL-RINSE **88**
OPERATING MINUTES TO END


9. Completed; display shows:

**OPEN
DOOR**

ON STANDBY
HANDS FREE START
BEDPAN & URINE BOTTLE
OPEN-DOOR
COMPLETED

ON STANDBY
HANDS FREE START
URINE BOTTLE BOWL
OPEN-DOOR
COMPLETED


To repeat or run an alternative cycle repeat steps 2, 3, and 4

 When opening door if the temperature in the chamber is above lockout temperature; display will show:




To close door without running cycle

After the door is opened; display shows:




- Manually press DOOR button
- Software will check that racks are in correct position before door closes. If racks have been moved, 'WRONG RACK' will flash. Place rack or racks into correct position for any of the 3 cycles and press DOOR button
- Door will close and go to standby


If racks are incorrectly positioned; display shows:



If racks are in correct positions; display shows:



After door is closed without running a cycle; display shows:



4.1 Release of Processed items

On completion of the cycle, it is recommended that all processed items are unloaded and visually inspected as clean as per the requirements AS 5369. It is recommended all items are removed on completion of the cycle and stored in a designated area to prevent any risk of contamination or mixing with unprocessed items.



If the processing status of a load or item is unknown or in question, it should be treated as unprocessed and not released.

5.0 Maintenance



All maintenance, preventative or breakdown shall be carried out by a qualified person. Failure to comply with this condition may result in unsafe conditions.

Preventative maintenance must be completed as per instructions below. Failure to comply may result in warranty claims being rejected.

The Malmet Bedpan/Urinal Bottle and Utensil/Bowl Washer Disinfector is self-cleaning; however proper care should be taken to ensure that the device is cleaned and maintained in accordance with maintenance instructions for Malmet Bedpan/Urinal Bottle and Utensil/Bowl Washer Disinfector, regulatory requirements and common-sense practices.

Preventative Maintenance Schedule

5.1 Daily Maintenance (Operator or Maintenance Technician)

- Run the Bedpan urine cycle to bring the device to operating temperature and disinfect the internal wash chamber area (do not load any items).
- Wipe out the inside of the door and gasket with warm water and detergent. A wipe with disinfectant is also desirable.
- Wipe over external stainless-steel panels with a stainless-steel cleaner (do not get cleaner into the opening of the USB port).
- Wipe the front control panel with a soft cloth and mild detergent as necessary. Care should be taken not to damage the digital display or to activate a cycle (turn off power at control panel).
- Check level in detergent bottle, replenish as necessary.
- Visually inspect for signs of leaking fluid from the device, as these may pose a pathogenic risk. Inspect around the chamber door, beneath the device and service connection points. If a spill or leak is observed, clean up as per the facility's infection control procedure and notify a maintenance technician to fix the cause of the leak.

5.2 Bi-Monthly Maintenance (Maintenance Technician)



WARNING 240/415 VOLTS!

ISOLATE DEVICE FROM ELECTRICAL SUPPLY BEFORE SERVICING



HOT SURFACES!

ALLOW DEVICE TO COOL PRIOR TO COMMENCING SERVICE WORKS

- Inspect for steam or fluid leaks, tighten unions, hose clamps and glands where necessary.
- Check flush nozzles for free rotation.
- Check spray nozzles are not blocked, insert a small pin into the orifice of the spray nozzle to ensure no blockages are present. If a build-up or blockage is detected have technician clean nozzles as per nozzle cleaning instructions below to ensure adequate performance.
- Remove level probe in steam generator tank and water tank, clean off any residue build up.
- Remove temperature probes in steam generator and chamber and clean off any residue build up.
- Visually inspect build-up of residue in steam tank, especially in locations with poor water quality.

- g) Check filter in the water inlet solenoid valve and clean as necessary.
- h) Check all electrical connections, and tighten if necessary.

Nozzle Cleaning Instructions

Remove and clean the sprays in sequence (do not mix sprays with other spray groups).

Chamber top to bottom: Cool down sprays x 2, Top sprays x 4, Bottom sprays x 4, Pan Tip sprays x 2 and rotary nozzles x 2.

Door: Urine bottle sprays x 4 – check that the holes are clear and clean as necessary.

Hold under tap and pressurise through the nozzle outlet in opposite direction of normal flow or clean in ultrasonic cleaner if available. An appropriately sized Oxy-Acetylene tip cleaner pin may be used to clean out spray nozzle orifices. In areas with hard water or high minerals in water supply, chemical de-scaler may be the quickest and easiest means to remove build-up of deposits. Replace before removing next spray group.

Stainless Steel Maintenance/Care

Under normal usage, stainless steel products require regular cleaning with a soft clean rag moistened with a mild detergent followed by a water-moistened clean rag and then a dry rag.

The #4 satin finish steel should be protected against Muriatic acid and caustic or abrasive materials and harsh cleaning detergents. In the event such agents cause discoloration, polish with a stainless-steel cleaner such as 3M Stainless Steel Cleaner & Polish and a pad.

5.3 Recommended Preventative Maintenance Schedule (to be performed by qualified maintenance personnel)

It is recommended that preventive maintenance is performed by a qualified maintenance technician every **3 Months** for devices in high use, **6 Months** for medium use or **Annually** for low use devices



WARNING 240/415 VOLTS!

ISOLATE DEVICE FROM ELECTRIC SUPPLY BEFORE SERVICING



HOT SURFACES!

ALLOW DEVICE TO COOL PRIOR TO COMMENCING SERVICE WORKS

1. Remove top and side panels and front bottom panel behind detergent door.

Note: Panel removal

- i) Remove 4 screws on the top panel
 - ii) Lift the side panel up to remove
 - iii) Remove 2 screws from bottom front panel
2. Remove and clean the sprays as per nozzle cleaning instructions on page 27.
 3. Check condition of rack positioning magnets.
Note: Replace magnets annually
 4. Check hinge screws for tightness on pan rack and bowl rack.
 5. Check condition of door gasket.
 6. Steam generator: Check for scale build up inside element housing and water tank. Remove and clean temperature probe and water level sensors.
 7. Water Tank: Check for scale build up inside tank and on element. Remove and clean temperature probe and water level sensor.

8. Remove and clean temperature probe at steam outlet.
9. Water inlet solenoids: Check filters for material build up, clean if necessary.
10. Make sure the water fill solenoids are completely shutting off and levels are not creeping up, if so, clean and/or replace the solenoid.
11. Detergent pump: Check hose on pump and hoses in detergent lines for cracks or signs of leaks.
12. Load Cell detergent bottle: Remove detergent bottle, unscrew load cell mounting plate. Ensure that the wire to the load cell is not strained when inspecting the load cell. Clean any detergent residue from under plate. Check load cell for corrosion, lightly spray with WD40.
13. Before replacing panels run purge reset cycle (see 5.18) three times to clear any airlocks in spray lines and check for leaks. Dust down components if necessary. Replace all panels and screws.
14. Cycle verification: Run the 3 cycle options starting with Bedpan/Urinal bottles. Download logged data to PC; and print a report of the last 3 cycles. Check temperatures attained and wash cycle times are within the operating parameters.
15. Replace battery on processor PCB every 12 months. Check time and date are correct, adjust if required.

Malmety will make available on request circuit diagrams, component parts lists, descriptions, calibration instructions, or information which will assist the user's appropriately qualified technical personnel to repair those parts of the device.

5.4 Post Maintenance or Repair Safety Checks

After service or repair; the device shall be checked to confirm correct cycle operation and that the device is in a safe operating state before being returned to service.

5.5 Door Safety Test

To check the correct functionality of the door safety and position sensors, open and close the door using the manual door open button. During operation the device will conduct a self-check for correct operation of these switches. Once complete, check door is firmly closed and is flush with outer panel and no faults are displayed.

5.6 Cycle Operation Check

After completion of the required checks above, a cycle operation check is to be completed. Run the device through a complete cycle ensuring correct sequence of operation as per 4.0 Cycle of Operation. Inspect for fluid leaks from around the chamber door, rear service connections and beneath the device. If no faults are found during the inspection, and the device displays no faults at end of the cycle, the device is fit for return to operation.

5.7 Post Maintenance and Repair of Electrical Systems

Additional to the above safety tests if any repair or service work has been carried out that could affect the electrical safety of the device; it must be inspected and tested as per the requirements of AS/NZS 3760 prior to being returned to service.

Information on replacement of these devices can be found in the Service Technicians Manual.

If device fails any of the above tests it is not safe for operation, and cannot be put into service until the cause of the fault has been rectified and successful completion of tests.

5.8 Replacement of Safety Devices

The following safety devices must be replaced and tested by a Malmety trained technician, failure to do so may impair the protection by the device.

- Door Actuator
- Door Micro Switch
- Water Tank Element Thermal Cut-out
- Steam Generator Thermal Cut-out

Information on replacement of these devices can be found in the Service Technicians Manual.



5.9 Validation

It is recommended that Validation of the device is carried out in accordance with the requirements of AS 5369 & ISO 15883, consisting of 3 identified stages.

IQ – Installation Qualification

OQ – Operational Qualification

PQ – Performance Qualification

Installation Qualification

Performed upon installation of equipment to ensure equipment has been installed in accordance with applicable Standards and manufacturer's instructions.

Operation Qualification

Performed immediately after installation, relocation, change of a service to the device or after repair prior to device being put back in service in accordance with applicable Standards.


Performance Qualification

Performed immediately after IQ or OQ for new installs or relocated equipment, or change of a service to the device or repair that might adversely impact the quality of the RMD.

5.10 Requalification

It is recommended requalification of the process is performed annually in accordance with the requirements of AS 5369 and ISO 15883 applicable parts.

5.11 Device faults

In the event of a fault or malfunction with the device during operation, abort the cycle by pressing the Power  On/Off button on the control LCD. Once the device is off, isolate the mains power supply and both the hot & cold-water services. Notify maintenance or service provider of the fault and ensure device is not used until the fault is rectified.


5.12 Fault Indication

The operation of the WDS is controlled by a micro-processor. The processor has fault detection capability and indicates faults by code on the LCD display as well as an audible buzzer.

In the event of a failure during cycle the device will go into fault. The display will show the selected cycle the device was operating in when the failure occurred along with a fault code corresponding to the cause of the fault. Access to the load will not be allowed until the cause of the fault has been fixed and a cycle has been successfully completed. Please refer to the Fault Codes section of this manual for a description of the fault code displayed.

FLT XXX flashing

Audible Buzzer will beep



After fault has been rectified a Purge/Reset will be required to reset the device (see 5.16). On completion of the “Purge Reset” device will restart previous cycle and will not allow access to the load until cycle has been successfully completed.

5.13 Fault Codes

FAULT CODE	STATE	DESCRIPTION	POSSIBLE CAUSES
901	System	Blockage	P or S trap blockage to drain
902	System	Door open during wash cycle	Faulty door micro-switch J11, J10, actuator, transformer fuse
903	System	Internal system communications error	Display unplugged, LINK J34 removed, loose processor PCB
904	System	Power interrupted during cycle	Mains power supply interrupted
905	System	Power button cycle aborted	Wash cycle aborted prior to completion
906	System	Blockage during purge reset	P or S trap blockage to drain
907/908	System	Steam generator temp sensor fault	Open/Closed Circuit
909/910	System	Tank water temp sensor fault	Open/Closed Circuit
911/912	System	Chamber temp sensor fault	Open/Closed Circuit
913/914	System	Chamber validation temp sensor fault	Open/Closed Circuit
915	System	Disinfection validation timer error	Processor Fault or Real Time Clock Error
101	Power On	Water fill timeout - Water tank and steam generator	Water level sensor failed or water supply low flow or turned off, water fill solenoids failed
102	Power On	Steam fill timeout	Steam fill solenoid failed, steam level sensor failed, blocked inlet hose
120	Power On	Blockage at power on	P or S Trap blockage to drain. Faulty sensor
141	Open Door	Door overcurrent trip	Door is obstructed, actuator motor faulty, door hinges seized
143	Open Door	Door timeout	Actuator motor faulty, blown transformer fuse, door hinges seized
145	Open Door	Chamber Temp Door Interlock	Chamber Temp above safe opening limit, faulty temperature probe, element stuck on.
149	Open Door	Door interlock switches fault	Door interlock switch J11 stuck on when door opens
181	Close Door	Door overcurrent trip	Door is obstructed, actuator motor faulty, door hinges seized



FAULT CODE	STATE	DESCRIPTION	POSSIBLE CAUSES
183	Close Door	Door timeout	Actuator motor faulty, blown transformer fuse, door hinges seized
185	Close Door	Chamber Temp Door Interlock	Chamber Temp above safe opening limit, faulty temperature probe, element stuck on.
186	Door Operation	Door Overcurrent Trip	Door Obstructed or Actuator Motor Fault
188	Door Operation	Door Timeout	Faulty Door Actuator or Micro Switch
190	Door Operation	Chamber Temp Door Interlock	Chamber Temp above safe opening limit, faulty temperature probe, element stuck on.
220	Cold Flush Fill	Water level sensor fail	Dirty conductivity sensor probe, faulty connections, faulty conductivity module
221	Cold Flush Fill	Water fill timeout	Water level sensor failed or water supply low flow, water fill solenoid failed, water supply off
241	Cold Flush Run	Low water level during wash	Check cycle time, pump rate too high, water Leakage
242	Cold Flush Run	Door open	Door opened during wash cycle, faulty door interlock switch J11 J10
243	Cold Flush Run	Pump run without solenoid	Check cycle times parameters, faulty solenoids
244	Cold Flush Run	Timeout before level reached	Low pump water flow, blocked rotary nozzles, faulty 1" solenoid, panel lockout switches J8 open circuit. Faulty level probe. Door interlock switch J11 J10
260	Cold Wash Fill	Water level sensor fail	Dirty conductivity sensor probe, faulty connections, faulty conductivity module
261	Cold Wash Fill	Water fill timeout	Water level sensor failed or water supply low flow, water fill solenoid failed, water supply off
281	Cold Wash Run	Low water level during wash	Faulty 4 -way solenoid, panel lockout switches J8 open circuit. Faulty level probe. Incorrect spray nozzles fitted
282	Cold Wash Run	Door open	Door opened during wash cycle, faulty door interlock switch J11 J10
283	Cold Wash Run	Pump run without solenoid	Check cycle times parameters, faulty solenoids
284	Cold Wash Run	Timeout before level reached	Low pump water flow, blocked spray nozzles, faulty 4 -way solenoid, panel lockout switches J8 open circuit. Faulty level probe. Door interlock switch J11 J10
300	Hot Wash+Det Fill	Water level sensor fail	Dirty conductivity sensor probe, faulty connections, faulty conductivity module
301	Hot Wash+Det Fill	Water fill timeout during detergent	Water level sensor failed or water supply low flow, water fill solenoid failed, water supply off
302	Hot Wash+Det Fill	Detergent timeout	Empty detergent bottle, blocked suction line/tubing, faulty flow sensor, faulty connections, faulty or blocked detergent pump
304	Hot Wash+Det Fill	Water heating timeout	Failed heater element, heater circuit breaker tripped, faulty temp probe. Element over temp cutout tripped, panel lockout switches J8 open circuit. Low water temperature, door interlock switch J11 J10
305	Hot Wash+Det Fill	Water heating low level	Water leakage, faulty wash solenoid, faulty level probe. Faulty connections, faulty conductivity module
321	Hot Wash+Det Run	Low water level during wash	Faulty 4 -way solenoid, panel lockout switches J8 open circuit. Faulty level probe. Incorrect spray nozzles fitted
322	Hot Wash+Det Run	Door open	Door opened during wash cycle, faulty door interlock switch J11 J10






FAULT CODE	STATE	DESCRIPTION	POSSIBLE CAUSES
323	Hot Wash+Det Run	Pump run without solenoid	Check cycle times parameters, faulty solenoids
324	Hot Wash+Det Run	Timeout before level reached	Low pump water flow, blocked spray nozzles, faulty 4 -way solenoid, panel lockout switches J8 open circuit. Faulty level probe. Door interlock switch J11 J10
340	Hot Rinse Fill	Water level sensor fail	Dirty conductivity sensor probe, faulty connections, faulty conductivity module
341	Hot Rinse Fill	Water fill timeout	Water level sensor failed or water supply low flow, water fill solenoid failed, water supply off
344	Hot Rinse Fill	Water heating timeout	Failed heater element, heater circuit breaker tripped, faulty temp probe. Element over temp cutout tripped, panel lockout switches J8 open circuit. Low water temperature, door interlock switch J11 J10
345	Hot Rinse Fill	Water heating low level	Water leakage, faulty wash solenoid, faulty level probe. Faulty connections, faulty conductivity module
361	Hot Rinse Run	Low water level during wash	Faulty 4 -way solenoid, panel lockout switches J8 open circuit. Faulty level probe. Incorrect spray nozzles fitted
362	Hot Rinse Run	Door open	Door opened during wash cycle, faulty door interlock switch J11 J10
363	Hot Rinse Run	Pump run without solenoid	Check cycle times parameters, faulty solenoids
364	Hot Rinse Run	Timeout before level reached	Low pump water flow, blocked spray nozzles, faulty 4 -way solenoid, panel lockout switches J8 open circuit. Faulty level probe, door interlock switch J11 J10
401	Steam Heatup	Steam water level timeout	Water leakage, steam fill solenoid failed, steam level sensor failed, blocked inlet hose
402	Steam Heatup	Steam Water Temp not maintained	Cool Down Spray or Solenoid failure/leak
403	Steam Heatup	Steam Water low level during heatup	Water leakage, steam fill solenoid failed, steam level sensor failed, blocked inlet hose
404	Steam Heatup	Steam water temperature high	Faulty temperature sensor, faulty heater relay, faulty contactors
405	Steam Heatup	Steam timeout	Faulty steam heater element, faulty steam heater relay, heater circuit breaker tripped, over temp switch tripped. Faulty temp probe, faulty contactors. Panel lockout switches J8 open circuit, door interlock switch J11 J10, leaking 4 – way solenoid valve.
421	Steam Disinfect	Disinfect water low level	Water leakage, , steam level probe failed
422	Steam Disinfect	Disinfect water temperature high	Faulty temperature sensor, faulty heater relay, Faulty contactors
423	Steam Disinfect	Disinfect timeout	Faulty steam heater element, faulty steam heater relay, heater circuit breaker tripped, over temp switch tripped. Faulty temp probe chamber steam outlet, faulty contactors. Blocked condenser coil. Blocked mesh on steam outlet, panel lockout switches J8 open circuit, door interlock switch J11 J10
424	Steam Disinfect	Temp Probe Differential	Differential between control and validation temp probes has exceeded tolerance.
425	Steam Disinfect	Disinfection	Difference between Control and Monitoring timebase > 3 second

FAULT CODE	STATE	DESCRIPTION	POSSIBLE CAUSES
801	Purge flush	Low water level during wash	Low pump water flow, blocked spray nozzles, faulty solenoid, panel lockout switches J8 open circuit. Faulty level probe. Incorrect spray nozzles fitted
802	Purge flush	Door open	Door opened when attempting pure rest, faulty door interlock switch J11 J10
803	Purge flush	Pump run without solenoid	Check cycle times parameters, faulty solenoids
804	Purge flush	Timeout before level reached	Low pump water flow, blocked spray nozzles, faulty solenoid, panel lockout switches J8 open circuit. Faulty level probe. Door interlock switch J11 J10
805	Purge Cold Fill	Water level sensor fail	Dirty conductivity sensor probe, faulty connections, faulty conductivity module
806	Purge Cold Fill	Water fill timeout	Water level sensor failed or water supply low flow, water fill solenoid failed, water supply off

5.15 Making equipment safe after incomplete operating cycle

If an incomplete cycle has occurred, it is assumed that the load being processed has not been subjected to a complete cycle as per the set process parameters. For this reason, access to the load is restricted and will not be allowed until a full cycle has been completed. At the end of an incomplete cycle or on powering up the device after an incomplete cycle, a fault code will display relating to the cause of the incomplete cycle (Please see section 5.13 Fault Codes for description of fault code). Depending on the nature of the fault a purge reset option may be available; this can be done by following the instructions as per section 5.16 of this manual. If the device does not allow a purge reset to be completed, or the fault re-occurs after a purge reset, it is recommended that the device is isolated and service is arranged to correct the fault.


5.16 Purge Reset (Operator)

1. Turn power OFF  at control LCD display
2. Turn power ON  at control LCD display
3. The display will alternate between the current fault/s and "PURGE", press SELECT  to start the Purge Reset cycle


When purge is complete the previous cycle will automatically start to ensure the load is disinfected, access to the load is restricted until a cycle has been completed to the required parameters

If the fault reoccurs contact maintenance or service provider

5.17 Purge Detergent Line (Operator)

1. Turn power OFF  at control LCD display
2. Turn power ON  at control LCD display

On powering back on the display will show TEST 302












3. The detergent line will automatically purge if necessary for a predetermined time or until the detergent is detected
4. Once detergent is detected the fault will clear and the previous cycle will be restarted

If the detergent fault is not cleared check detergent level and repeat the purge process again, if the detergent fault continues and is unable to be cleared contact maintenance or service provider

The following functions must only be undertaken by Malmet service personnel or Malmet-trained facility maintenance personnel.

5.18 Manual Purge Reset (Maintenance)

1. Turn power OFF  at control LCD display
2. Turn power ON  at control LCD display
3. Hold down both SCROLL  and SELECT  buttons, when MENU is displayed, release buttons
4. SCROLL  to SERVICE menu and press SELECT 
5. SCROLL  to RESET menu and press SELECT  'CONFIRM' should then display
6. Press SELECT  again to start Purge Reset cycle

When purge reset is completed the display will resume normal operation













5.19 Manual Purge Detergent Line (Maintenance)

1. Turn power OFF  at control LCD display
2. Turn power ON  at control LCD display
3. Hold down both SCROLL  and SELECT  buttons, when MENU is displayed, release buttons
4. SCROLL  to TEST menu and press SELECT  RELAY
5. SCROLL  to DETERG menu and press SELECT  START press SELECT  START/RUNNING

The detergent pump will run for 5 seconds, if detergent line completely purged display shows PASS; then DETERG

If display shows FAIL, check detergent level and repeat Purging of the Detergent Line







5.20 Safe Access to the Load or Chamber in Fault Condition

1. Turn power OFF  at control LCD display
2. Turn power ON  at control LCD display
3. Hold down both SCROLL  and SELECT  buttons, when MENU is displayed, release buttons
4. SCROLL  to TEST menu and press SELECT  RELAY
5. SCROLL  to DOOR OPEN and press SELECT  to turn on the door opening relay (ensure hands are clear and no obstructs or objects are in the opening path of the door)
6. Once the door is open, press SELECT  to turn off the door open relay
7. To close the door, press SCROLL  to DOOR CLOSE and SELECT  to turn on the door closing relay (ensure hands are clear and no obstructs or objects are in the opening path of the door)
6. Once the door is closed, press SELECT  to turn off the DOOR CLOSE relay



Treat all items in the load/chamber as unprocessed, do not use until they can be cleaned or processed accordingly.

To EXIT Service Menu:

1. SCROLL  to EXIT press SELECT 
2. SCROLL  to EXIT press SELECT 
3. SCROLL  to EXIT press SELECT 

Once all menus have been exited the LCD screen will turn blank, once unit is powered on from this state it will return to normal operation.



6.0 Technical Data

6.1 Power and Water Consumption

Power and Water Consumption Data (for all WDS 1ph models)						
Type of Cycle	Disinfecting Time/ Temp	Avg Cycles per/Hr	Avg Cycle min/sec	Avg kWh	CW Avg Lt Per cycle	HW Avg Lt per cycle
Bedpan / Urinal Bottle	1 min @ 90°C	4.7	12.45	0.62	35.6	21.6
Urinal Bottle Only	1 min @ 90°C	5.3	11.15	0.572	24.4	20.7
Bowl / Utensils	1 min @ 90°C	5.2	11.30	0.572	24.4	20.7

Power and Water Consumption Data (for all WDS 3ph models)						
Type of Cycle	Disinfecting Time/ Temp	Avg Cycles per/Hr	Avg Cycle min/sec	Avg kWh	CW Avg Lt Per cycle	HW Avg Lt per cycle
Bedpan / Urinal Bottle	1 min @ 90°C	5.1	11.45	0.60	35.6	21.6
Urinal Bottle Only	1 min @ 90°C	5.8	10.15	0.574	24.4	20.7
Bowl / Utensils	1 min @ 90°C	5.7	10.30	0.574	24.4	20.7




NOTE: Values may change due to operating and supply service conditions




6.2 Device Specifications

Model WDS Single Phase	Electrical Rating	Volts	240V	APPROVALS
		Phase / Hz	1 ph / 50 Hz	
		Amps	20 Amps	
	Heating Elements	Steam Generator	4500W 240V 3 x 6.25A Star	
		Water tank	4500W 240V 3 x 6.25A Star	
Model WDS Three Phase	Electrical Rating	Volts	415V	
		Phase / Hz	3 ph / 50 Hz	
		Amps	20 Amps	
	Heating Elements	Steam Generator	4500W 240V 3 x 6.25A Star	
		Water tank	7200W 240V 3 x 10A Star	
	Heating Element Contactors	Steam heating element	3 pole 7.5kw	
		W/Tank heating element	3 pole 7.5kw	
ELECTRICAL (Common components)				
Element over temperature cut-out		Steam Generator	Capillary bulb thermal cut out Manual reset 115°c 25A	
		Water tank	Capillary bulb thermal cut out Manual reset 100°c 25A	
Fuse: Transformer		Cylinder type 5 x 20mm	T 5A H 250Vac Time Delay	
Door Actuator		Type LA12 Self-Locking	24VDC	
Interlock Micro Switch		SPDT	240VAC 16A	
Processor module battery		Coin Cell	CR2032 Lithium 3Volts	
PCB HARDWARE				
PCB (Printed Circuit Boards)		Display LCD	D134411	
		Heating switch module	D132213	
		Control module	D134412	
		Processor module	D133447	
		Conductivity sensor	D134418	
PC Interface Port		USB Type B		
Environment conditions (Operating)		Location	Indoor Use	
		Altitude	<2000m	
		Temperature	+10°C to +25°C	
		Relative Humidity	+30% to 75%	
		Mains voltage fluctuations	240±10%	
		Overvoltage	Category II	
		Pollution	Degree II	
WASH CHAMBER				
Loading Configurations		Bedpan Cycle	4 x Urinal Bottles and <u>One (1)</u> of the following:	
			2 x Small slipper pans & Lid	
			2 x Large Slipper pans & Lid	
			2 x Standard Bedpans & Lid	
		Urinal Cycle	4 x Urinal bottles Only	
		Bowl / Utensil Cycle	2 x Large bowls	
			3 x Medium bowls	
			6 x Large kidney dishes	
			6 x Medium kidney dishes	
			Assortment of Smaller bowls and kidney dishes using optional accessory basket.	
Chamber Specifications		Material	1.2mm 304 Stainless Steel	
		Design Life	Minimum of 10,000 Cycles	
		Usable Volume	100L	
		Total Volume	115L	
WASH SYSTEM				
Nozzles and Sprays		Flush nozzles	2 x Rotary Nylon	
		Tip sprays	2 x WL1/2 120° 1/8 BSP SS	
		Top sprays	4 x WL1/4 90° 1/8 BSP SS	
		Bottom sprays	4 x WL1/2 60° 1/8 BSP SS	



	Urine Bottle sprays	4 x WL1/4 90° 1/8 BSP SS		
	Cool Down sprays	2 x LW1 120° 1/4 BSP SS		
Wash pump	Horizontal multistage	240V 1ph 3.1A		
Detergent	Malmet specific (5Lt)	Detergent Concentrate Caustic Alkaline	ARTG Class 1	
Detergent pump	Self-priming peristaltic pump	240V 50Hz		
PLUMBING				
Soil line Connection	Type	‘S’ or ‘P’ Trap (6mm PE) 80°	Complies to AS 2887-1993 (R2017) AS/NZS 3500.1:2018 & AS/NZS 3500.2:2018	
	Size	100mm ID		
Backflow Prevention	Air Gap (RAG)	Water storage tank	AS 2845.2 - 2010	
Cold water inlet Supply Potable Water (Final Rinse Water)	pH	6.5-8.5pH		
	Chloride	<30mg/L		
	Water Hardness	<2.5mmol/L		
	Water conductivity	<850 µS/cm		
	Temperature	15 – 25°C		
	Pressure	100 – 350kpa		
	Nominal Flow Rate	10L/Min		
	Solenoid Valve 1 in 2 out	240V - GB¾ male	WRAS Certified	
Hot Water inlet supply	Temperature	55 – 60°C		
	Pressure	100 – 350kpa		
	Nominal Flow Rate	10L/Min		
	Solenoid Valve 1 in 1 out	240V - GB¾ male	WRAS Certified	
Hose – valve to Water Tank	S/Steel Braided	10mm	WRAS Certified	
Hoses – Mains Water Supply to Inlet Valve	S/Steel Braided	G¾ Hex Nut Each End 1.5M x 10mm	Certified to AS/NZS 3499-2006 (R2016)	
CONSTRUCTION MATERIALS				
Materials	Frame	25 x 25 x 1.6 tube SS		
	Wash chamber	304 SS grade 4		
	Door	304 SS grade 4		
	External panels	304 SS grade 4		
	Water tank	316 SS grade 2B		
	Steam generator	316 SS grade 2B		
	Fasteners	304/316 SS		
	Thermal insulation chamber and water tank	10mm Cello HR 290 VLC HT		
Principal Heavy Component	Main Wash Pump	11kg		
TRANSPORTATION				
HANDLING & STORAGE CONDITIONS	 Fragile	 Keep away from rain	 Do not stack	
	Temperature	-5°C to +50°C		
DIMENSIONS & WEIGHTS				
Weights	Nett:	Operating:	Shipping:	Shipping (crated):
	140kg	176kg	165kg	217kg
Dimensions (H x W x D)	1695X610 x 625		1870 x 810 x 800	1910 x 810 x 800
Floor Loading at Each Support as viewed from front of unit. * Loadings are a guide only & will vary between installations.	Front Left	Front Right	Rear Left	Rear Right
	58kg*	44kg*	49kg*	25kg*

6.3 Wiring Diagram (1ph devices) Manufactured before June 2020





Warranty Statement

This warranty is provided, and operates in addition to, the statutory warranties Malmet (Australia) Pty Ltd ("Malmet") provides to any consumer under the Australian Consumer Law (if applicable) or by virtue of any other applicable legislation.

Subject to the following conditions, we provide, from the date of purchase, the following warranty on Malmet devices and spare parts for products manufactured by Malmet and sold in Australia:

- Functional components found within the device to be defective in workmanship or material will be repaired or replaced free of charge subject to the periods of warranty specified in the table below.
- A decision regarding whether the defective components will be repaired or replaced will be determined at the sole discretion of Malmet or its authorised agents or representatives.
- The structural warranty covers any structural components within the device, which fail to perform their intended function due to faulty manufacture or deterioration within the warranty period.
- Parts replaced in devices under warranty are warranted for the balance of the original warranty period for that device.

Malmet Devices	
Device Components	Parts & Labour
Structural Guarantee	2 Years from Date of Purchase
All other components	2 Years from Date of Purchase

Malmet Spare Parts	
1 Year from Date of Purchase	

The installer is responsible for the correct installation, start up and demonstrating the operation of the product. They are also responsible for issuing the relevant certificates of compliance (these may differ from state to state).

CONDITIONS AND EXCLUSIONS

- Device must be installed and commissioned according to Malmet's instructions (outlined in Malmet Operation, Maintenance and Installation Manual) and operated to the purpose it was designed.
- Device must be serviced as instructed in the Operation, Maintenance and Installation Manuals.
- To the extent permitted by law, this warranty shall not cover damage, malfunction or failure resulting from accident, misuse or misapplication, improper or unauthorised repair, neglect or modification or use of unauthorised replacement parts or accessories, inclusive of detergent, or improper voltage. The warranty may be void if the serial number is removed or altered.
- Parts damaged in transit back to Malmet Leeton due to poor packaging could result in warranty claim being rejected in part or in full.
- Any part tampered with or which has been altered by unauthorised repairs and/or modifications will be rejected under a warranty claim to the extent permitted by law (to the extent the Australian Consumer Law applies, Malmet will assess the extent to which the tampering or unauthorised repairs contributed to the failure).
- Reasonable access must be allowed for maintenance. If any additional equipment is needed to provide access to the device, this must be provided (and paid for) by the owner.
- It is the owner's responsibility to provide safe access to the device. Malmet, or any of its authorised service agents, may refuse to perform maintenance or warranty work if access is unsafe, as determined by Malmet or any of its authorised service agents acting reasonably.
- Should a warranty claim be rejected you will be advised in writing with a full explanation of our reasons.
- Malmet have a Warranty Claim Procedure that is fair to our customers and provides an efficient system of replacement and/or repair of faulty parts. If at any time you believe we are not meeting our commitment to you please contact Malmet Head Office via email: info@malmet.com.au

- To the extent permitted by law, no responsibility will be accepted for outside elements including, but not limited to storms, pest and vermin that may cause damage to the device.
- To the extent permitted by law, no responsibility will be accepted for damage incurred as a result of, or incidental to, electrical surges or brown outs or for any other consequential damages.
- If there is no certificate of compliance for plumbing or electrical, Malmet reserves the right to refuse service on non-compliant installations.
- To the extent permitted by law, claims for damage to contents, carpet, ceilings, foundations or any other consequential loss either direct or indirect resulting from, power spikes, incorrect operation, incorrect installation, faulty product or any other cause, are excluded.
- This warranty, and to the extent permitted by law, any warranties owed by Malmet under the Australian Consumer Law or other applicable legislation, are not transferrable and cannot be sold, assigned or transferred in any other way from the purchaser to any other person.
- To the extent permitted by law, unauthorised use of any parts that were not supplied or approved for use in the applicable device by Malmet will result in this warranty and any warranty claims applicable to that device being void.
- Warranty labour (service work) shall not include devices located outside of city metropolitan areas of Melbourne, Sydney, Adelaide, Perth and Brisbane. Costs outside these areas shall be borne by the owner. The owner shall be notified of this prior to the warranty call out.
- Warranty labour (service work) shall be performed during normal business hours (Monday – Friday 7am – 4pm), excluding public holidays.
- Warranty labour (service work) performed outside of normal business hours, shall be charged at Malmet's or its authorised representative or agent's standard after-hour labour rates.
- Warranty relating to spare parts covers parts only and does not include any associated labour costs.

To the extent permitted by law, a charge will be made for work done or a service call made where:

- There is no fault apparent with the device, as determined by Malmet or its authorised representative or agent acting reasonably.
- The defective operation of the device is due to failure of electricity or water supply.
- Defects are caused by neglect, incorrect application, abuse or by accidental damage of the device.
- An unauthorised person has attempted to repair the device.
- Harsh environmental situations including, but not limited to, water quality that may cause the water tank damage cannot be covered under this warranty



HOW TO MAKE A CLAIM UNDER THIS WARRANTY

If you believe there is a defect in a device you have purchased from Malmet, you must notify Malmet in writing of such defect, by sending an email (Notice of Defect) to info@malmet.com.au prior to the expiration of the applicable warranty period set out in this warranty.

For the avoidance of doubt, Malmet must receive your Notice of Defect prior to the expiration of the warranty period.

To the extent permitted by law, Malmet will not reimburse you for any expense you incur in claiming or attempting to make a claim for repair or replacement of a component under this warranty.

Please complete details below:

Date Purchased:	Warranty Expiry Date:
Sold To:	For Service Contact:

PROOF OF PURCHASE

Please retain your proof of purchase (receipt, invoice or commissioning certificate is accepted).

E.&O.E.

In the interest of continued product improvement, Malmet reserves the right to alter specifications without notice.

AUSTRALIAN CONSUMER LAW DISCLAIMER (APPLIES ONLY TO THE EXTENT YOU ARE A 'CONSUMER' WITHIN THE MEANING OF THE AUSTRALIAN CONSUMER LAW):

Malmet goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

Queensland Distributor

EVOCARE AUSTRALIA PTY LIMITED

A.B.N. 98 078 566 604

Trading as **EVOCARE** and **L&M EQUIPMENT**

P.O. Box 145, Everton Park Qld. 4053

Ph: 07 3355 8000

Website: <http://www.evocare.com.au>

Email: sales@evocare.com.au

workshop@evocare.com.au

warehouse@evocare.com.au

accounts@evocare.com.au



Manufactured by

Malmet (Australia) Pty Ltd

ABN 95 001 717 791

www.malmet.com.au