

Adaptable Safety Sides



Clinic | Care | **Homecare** | Assist



Bed insert Lippe IV with adaptable safety sides in the Relax senior bed

The adaptable safety sides from Burmeier give people in need of care the necessary protection they need from falls, without seriously restricting their freedom. Burmeier offers the adaptable safety sides for the **Arminia III 24 Volt, Dali II 24, Economic II/III** care beds and the bed insert **Lippe IV**.

The safety sides support the **mobility** of the occupant and enhances their **independence**. With a length of 120 cm it offers a higher level of safety, but at the same time leaves sufficient room to get out of bed.

In the higher position the side rails can be used as **handrails**, which provide for a more secure grip and also support mobilisation. The design is **open** and **clearly** structured. In the lowered position the safety sides are located underneath the mattress base ensuring that it is **more comfortable** to get in and out of bed.

Due to the easy operation, carers or relatives can operate the adaptable safety sides **without any problems** or **physical effort**.

Advantages at a glance

- Easy and fast assembly – even for retrofits
- Intuitive and easy to operate
- Can be combined with full-length safety sides on a care bed
- Protective height: 39 cm, length: 120 cm

 **Evocare**[®]
www.evocare.com.au
Ph: 07 3355 8000 Fax: 07 3355 5043
Email: sales@evocare.com.au



Arminia III 24 Volt care bed with a combination of adaptable and full-length safety sides

Adaptable safety sides for our Arminia III 24 Volt, Dali II 24 Volt and Economic II/III care beds and the Lippe IV bed insert

The adaptable safety sides can easily be fitted to the mattress base instead of full-length safety sides – even for retrofits. They are attached to the mattress base frame with 2 mounting brackets. The swivelling adaptable safety side is easily and continuously moved past the frame.

The adaptable safety side meets all the standard requirements for safety sides. At the same time it leaves enough room at the higher position for getting in and out. It follows the Burmeier philosophy of avoiding the restriction of freedom as much as possible.

Variant

right*

left*

*viewed from the foot end of the bed



Adaptable safety side in higher position (right)



Attached with mounting brackets



Adaptable safety side in lower position (right)

Burmeier GmbH & Co. KG
 Pivitsheider Straße 270, 32791 Lage
 Phone: +49 (0) 5232 9841 - 0
 Fax: +49 (0) 5232 9841 - 41
 Email: info@burmeier.com
www.burmeier.com

