ML125 SCHEDULE 8 DRUG REFRIGERATOR



FEATURES

- Complies with all the regulations of storing schedule 8 drugs.
- Bolts to the floor and back wall.
- Magnetic keypad lock with battery backup, will work even when the power goes out.
- Keypad lock can be used with 4-digit pin code or swipe card so you can limit who has access to the products.
- Solid door this means what products are in there cannot be seen
- High / Low Readout
- Door Ajar Alarm
- Battery back-up for controller
- Data Logger included
- Door Fan Switch
- Port Entry
- Fan forced refrigeration system, Automatic water disposal











ML125 SCHEDULE 8 **DRUG REFRIGERATOR**



MODEL	ML125P
Gross Capacity	125 L
External Dimensions in mm (w/d/h including castors or legs, width includes keypad)	660 x 610 x 840
Internal Dimensions in mm (w/d/h)	505 x 442 x 625
Packed Dimensions in mm (w/d/h including pallet)	760 x 705 x 910
Gross weight – kg (packed weight)	68
Net weight – kg	64
Voltage	240V
Temperature Range	2°C to 8°C
Ambient Temperature Range	10°C to 32°C
Refrigerant	R600
Cooling System / Defrost	Fan Forced
Exterior Cabinet Finish	Enamelled Steel
Inner Liner Material	ABS
Door	Solid
Door Hinge	RH
Door Handle	Pull Handle
Type of Controller	Digital
Temperature Display	External
Data Logger	Included
Door Ajar Alarms	Yes
High and Low Readout	Visual / Audible
Remote Alarm Contacts (BMS Connectivity for Alarms only) (BMS Connectivity for Temperatures)	Yes Optional
Port Entry	Yes
Battery Back-Up for Controller	Yes
Power Requirement	10amp GPO
Castors	No, Adjustable Levelling Legs
Light	LED
Noise – dB(A)	45





EVOCARE AUSTRALIA PTY LIMITED

A.B.N. 98 078 566 604

Trading as **EVOCARE** and **L&M EQUIPMENT**

P.O. Box 145, Everton Park Qld. 4053

Ph: 07 3355 8000



Website: http://www.evocare.com.au sales@evocare.com.au Email: workshop@evocare.com.au warehouse@evocare.com.au

accounts@evocare.com.au

We reserves the right to alter the specifications without prior notice.