

Clinic | Care | Homecare | Assist



The
paediatric cot
for special requirements





Safe and dependable

The Junior paediatric cot is distinguished by its exceptionally safe and dependable functions. Its practical features greatly facilitate the daily care routines performed by nursing staff.

Its sturdy, high-quality design makes the cot a safe place to sleep and at the same time a perfect place to play during the daytime. The cheerful colours enhance the children's acceptance of the cot and also promote visual perception.

The cot is suitable for small children of up to 100 cm in height and is designed for a safe working load of 50 kg.

With a design that meets all the relevant standards, the cot will provide many years of particularly stable and reliable service.



Ergonomic nursing

The height of the mattress base can be adjusted from 75 to 91 cm. This guarantees that an ergonomic height can be found to suit both hospital staff and family members. An easy release mechanism allows the cot to be adjusted to several intermediate heights (image 1 and 2).

Besides this height adjustment facility, the cot's mattress base can, also be tilted by hand to a 7° Trendelenburg and reverse-Trendelenburg position (large figure). This makes it the ideal cot for administering medical and nursing care to young patients.



Safety and protection

The easy to use safety sides of the Junior cot provide quick access to the child (image 2). At a height of 81 cm, they are also sufficiently high to protect children who can already stand up by themselves (image 3). It is also possible to adjust the safety sides to several intermediate positions in order to meet all nursing re-

quirements. The safety sides' release mechanism is designed to be operated from the outside only (image 1). This prevents children from lowering the safety sides themselves and potentially falling out of bed.



Accessories: practical, clean, flexible

The bed can be complemented with useful accessories. These optional extras are attached in a way that makes them quick and easy to reach. They include a pull-out laundry basket, a nappy bin and two sleeves for holding infusion stands. The laundry basket makes it easier for staff and family members to perform tasks such as changing the child's clothes or making the

bed since the basket can be used to hold or store clean or used laundry (image 2). The nappy bin enables used nappies to be disposed of quickly and hygienically (image 1). The sleeves for holding infusion stands also ensure that young patients can be provided with medication (large figure).

Technical data

	Junior
Dimensions and weights	
External dimensions	80 x 167 cm
Mattress base (mattress dimensions)	69 x 135 cm
Mattress base height	75 – 91 cm
Height of bed	170.1 cm
Safe working load	50 kg
Floor clearance of chassis	66 cm (safety sides at highest level)
	16 cm (safety sides at lowest level)
Total weight	77 + 5 = 82 kg (bed + accessories)
Protection height of safety sides	up to approx. 81 cm
Adjustment ranges	
Tilting to Trendelenburg position	7°
Tilting to reverse-Trendelenburg position	7°
Mattress base height	75 – 91 cm
	Height and Trendelenburg position adjustable with a locking mechanism
Materials	
Chassis and mattress frame	Welded steel structure with durable polyester powder coating
Accessory equipment and optional features	
	Pull-out laundry basket ○
	Nappy bin ○
	Location sleeve for infusion stand ○

○ = optional

Base colours

head and footboards, long sides, handles, infusion sleeves

Crystal white NCS 0500-N	Argentum plus
Traffic yellow RAL 1023	Sky blue RAL 5015

Accent colours

safety sides, mattress frame, lever for height adjustment, guide rail for laundry basket, backrests

Crystal white NCS 0500-N	Argentum plus
Traffic yellow RAL 1023	

Ihr Bett – unsere Leidenschaft
Your bed – our passion

EVOCARE AUSTRALIA PTY LIMITED

A.B.N. 98 078 566 604

Trading as **EVOCARE** and **L&M EQUIPMENT**

P.O. Box 1144, Stafford Qld. 4053

Ph: 07 3355 8000 Fax: 07 3355 5043

Website: <http://www.evocare.com.au>

Email: sales@evocare.com.au

workshop@evocare.com.au

warehouse@evocare.com.au

accounts@evocare.com.au



Technische Änderungen vorbehalten. Leichte Farbabweichungen der Originale sind drucktechnisch bedingt. Stand 03/2017. Version 2A. Artikel-Nr. 600557.
Subject to technical changes without notice. The shown product colours may differ slightly from the actual product. Last updated 03/2017. Version 2A. Item no. 600557.

Finden Sie Ihren Händler auf unserer Webseite:
Find your local partner at our website:
stieglmeyer.com/partner

Stieglmeyer GmbH & Co. KG

Ackerstraße 42, 32051 Herford, Germany

Phone +49 (0) 5221 185 - 0

Fax +49 (0) 5221 185 - 252

E-Mail info@stieglmeyer.com

www.stieglmeyer.com

