

Etac 2Move

Transfer board for horizontal transfer



Etac 2Move transfer board is used to bridge the gap and even out the height difference between two horizontal positions. It is used for horizontal transfer between beds, bed to shower trolley, bed to x-ray table, stretcher to bed, etc. It is simple to transport and store.

2Move transfer board is available in two lengths; short (90 cm) or long (180 cm). The long model is available in two variants; narrow foldable or narrow non-foldable.

2Move is always used in combination with a low friction cover that by sliding around the board, minimizes the friction and facilitates the transfer. Covers are available in two alternatives; waterproof nylon or disposable.



Disposable cover in roll.



Etac 2Move transfer board, including nylon cover.

Intended Use

- Bridge the gap when transferring between two horizontal positions
- Transfer between beds
- Transfer between bed and shower trolley
- Transfer between bed and stretcher, x-ray table, operating table

Description

Description	Width	Length	Item no.
2Move Transfer Board, Short	50 cm 19.7"	90 cm 35.4"	IM8050090
2Move Transfer Board, Long foldable	50 cm 19.7"	180 cm 70.9"	IM8050180
2Move Transfer Board, Long non-foldable	50 cm 19.7"	180 cm 70.9"	IM8150180

Covers

2Move Nylon Cover, waterproof, for IM8050090	50 cm 19.7"	90 cm 35.4"	IM4045
2Move Disposable Cover for IM8050090, roll with 150 pcs	60 cm 23.6"	100 cm 39.4"	IM4031
2Move Nylon Cover, waterproof, for IM8050180 & IM8150180	50 cm 19.7"	180 cm 70.9"	IM4044
2Move Disposable Cover for IM8050180 & IM8150180, roll with 75 pcs	60 cm 23.6"	200 cm 78.7"	IM4030

Specifications - Board

Top: Nylon Inside: Plywood
Handles: Polyester Cover: Nylon

Washing Instructions

Can be wiped off with disinfectant.



Specifications - Cover

Nylon

max 150 kg max 330 lbs

Washing Instructions

Wash as instructions.
Do not use fabric softener or bleach

