


Etac Transfer Boards

Mind the gap...

Etac Butterfly

Etac Butterfly is a long, wide and stable transfer board. It has cut outs for positioning against the wheel of the wheelchair and non-slip pads on the underside, which allows for good positioning and help prevent sliding out of place. Etac Butterfly has a larger sitting and sliding area, which is helpful when the transfer needs to be made in short steps.

Etac Butterfly can be used with the special Butterfly GlideWing, which is mounted on the board. The glide wing slides along the board and provides additional slide function during the entire transfer.



max 250 kg max 550 lbs

Etac Dolphin

Etac Dolphin is a sturdy and stable transfer board uniquely designed for transfer to the toilet. The board has cut outs for positioning against the wheel of the wheelchair and non-slip pads on the underside, which allows good positioning and avoids sliding out of place. The cut out allows the user to remain on the board while on the toilet.



max 150 kg max 330 lbs

Intended Use

- Transfer from bed to wheelchair / shower chair
- Transfer in and out of the car
- Transfer between wheelchair, chair, toilet, shower chair, commode

Description	Width	Length	Item no.
Butterfly, transfer board	32 cm 12.6"	66 cm 26"	IM401 ✓
Butterfly, glide wing L	40 cm 15.7"	31 cm 12.2"	IM4001 ✓
Butterfly, glide wing XL	52 cm 20.5"	31 cm 12.2"	IM4002 ✓
Dolphin Toilet Board, white		69 cm 27.2"	IM409 ✓

Specifications - Boards

Fiberglass

Washing Instructions

Transfer boards can be cleaned with soapy water or disinfectant.

Specifications - Glide Wings

Material: Nylon

Padding: Polyester

Washing Instructions

Wipe off with soapy water or disinfectant, washed as described below.

Do not use fabric softener or bleach.

