

Etac HandyGrip

Adjustable grip to aid turning in bed

Etac HandyGrip provides help with independent turning in bed. It can also be used when self transferring from wheelchair to bed. Soft, compact and easy to hide away. The adjustable buckle allows the length and position of HandyGrip to be varied, providing the user with a good grip at the required position.

Soft, compact and easy to hide away.



Intended Use

- Support independent turning in bed
- Support independent transfer to bed

Description

HandyGrip

Width

5 cm 2"

Length

160 cm 63"

Item no.

IM461

Specifications

Polyester

Belt Buckle: Acetal and nylon



Washing Instructions

Do not use fabric softener or bleach.

