

# Etac MultiGlide Glove

Etac MultiGlide Glove can be used in many different manual handling situations. Using the glove facilitates placing the hands between the user and mattress or cushion. Small re-positioning adjustments, movements and smoothing out the bed sheets are common activities for which the gloves are very useful, and will provide an improved solution for the user and the caregiver. The gloves are also useful for instance when checking the pressure of an air cushion on a wheelchair, whilst the user is still sitting on it.



**Intended Use**

- Positioning
- Small re-positioning adjustments
- Hand check of anti-decubitus cushions and mattresses ("Bottoming-Out" test)



Description	Width	Length	Item no.
MultiGlide Glove (1 pair)	20 cm 7.9"	50 cm 19.7"	IM35 ✓

**Specifications**  
Nylon with polyurethan coating.

**Washing Instructions**  
Do not use fabric softener or bleach.

