

# Etac OneWayGlide

OneWayGlide is a cushion that allows sliding easily in one direction whilst resisting sliding in the other direction. It is suitable for users that slide forward when sitting. OneWayGlide facilitates independent or assisted positioning back in the chair whilst helping to prevent the user from sliding back forward again.

## Etac OneWayGlide – Standard

Available in 3 variants; OneWayGlide nonslip, OneWayGlide velour, OneWayGlide pad. OneWayGlide non-slip is a tubular cushion with a non-slip outer material that grips to any surface. It can be combined with a cotton cover when necessary.

OneWayGlide velour is a tubular cushion with 'velour material. It is more comfortable to sit on, but needs to be used on a textile surface. OneWayGlide pad is a one layer pad with velour top suitable for use directly on a textile surface.



## Etac OneWayGlide – Long

Etac OneWayGlide Long is a long and open cushion with handles and velour top surface. The length allows for re-positioning to be performed several times in one application. The handles provide a good grip, allowing the caregiver to assist.



### Intended Use

- Positioning back in chair / wheelchair / day chair

Description	Width	Length	Item no.
OneWayGlide, tube non-slip S	37 cm 14.6"	43 cm 16.9"	IM39/8 ✓
OneWayGlide, tube non-slip M	45 cm 17.7"	45 cm 17.7"	IM38/8 ✓
OneWayGlide, tube non-slip L	80 cm 31.5"	60 cm 23.6"	IM37/8 ✓
OneWayGlide, tube velour	37 cm 14.6"	43 cm 16.9"	IM45/8
OneWayGlide, tube, velour	40 cm 15.7"	50 cm 19.7"	IM50/8
OneWayGlide, pad, velour	43 cm 16.9"	37 cm 14.6"	IM44/8
OneWayGlide Long	40 cm 15.7"	130 cm 51.2"	IM34/8 ✓

### Specifications

Padding: Polyester  
Inside: Nylon  
Material: Polyester & Anti-slip  
Material (Velour): Polyester, Nylon

### Washing Instructions

Do not use fabric softener or bleach.

