

Torkel – Toilet paper tongs

Torkel toilet paper tongs provide a functional grip and extended reach to handle intimate hygiene.

Easy to handle

When grip and mobility function are restricted, these tongs can be used for toilet paper, wash cloths and small towels. They are made of a material that can withstand boiling.

How to use Torkel toilet paper tongs:



1. Tear off a piece of toilet paper (approx. 40 cm). Open the tongs by squeezing the handle, then thread the paper between the tongs.



2. Roll the paper loosely around the tongs so that it covers the tip.



3. Place the tongs between the legs and wipe.



4. Open the tongs and shake to drop the used toilet paper into the toilet. Rinse the tongs.



Description

Torkel toilet paper tongs

Length

33 cm 13"

Weight

100 g 3.5 oz

Material

Polyamide plastic

Item no.

8030 4001



max 85°C max 185°F
for maximum of 3 minutes



Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 100°C / 212°F (boiling) if so required.