

Socky – Stocking aids

Socky stocking aids make it possible for people with impaired mobility to put on socks, stockings and support stockings. This range is particularly suitable for users with back, hip or knee problems.



Socky stocking aid, short

Socky short stocking aid is small, flexible and easy to use. The narrow front in combination with the nylon fabric makes it is easy to pull the sock onto the stocking aid. The foot slides easily on the fabric bringing the sock into place.


Socky stocking aid, long

Socky long stocking aid has the same features as the short version, and includes a handle and plastic sheet on the underside. The sheet slides easily on carpets and all floors while the grip makes it easier to pull up the stocking. For tights, use two Socky stocking aids connected with the linking device.

Material

Fabric: Polyamide

Plastic parts: Polypropylene



Socky stocking aid, short


Socky stocking aid, long

Description	Length	Weight	Item no.	
Socky stocking aid, short	40 cm 15¾"	94 g 3.3	8060 1002	✓
Socky stocking aid, long	60 cm 23½"	146 g 5.1 oz	8060 1001	✓
Linking device			8060 1021	✓



Socky stocking aid for support and compression stockings

This stocking aid is designed for putting on support stockings of various lengths and types. It has large handles which are also helpful for an assisting carer.



Description	Length	Weight	Item no.	Material
Socky stocking aid for support and compression stockings	100 cm 39¾"	90 g 3.2 oz	8060 1003	Fabric: Polyamide Plastic parts: Polypropylene