

# Etac SmartMove

## Support pad with handles

Etac SmartMove acts as extension of the arms, allowing the caregiver to offer support and assistance whilst maintaining a good ergonomic position and grip.

Etac SmartMove pad can be used to support and help transfer from sitting to standing, standing to sitting, to move higher up in bed or move back / forward in a chair. SmartMove provides the caregiver with a place to grip and to maintain a good working position. The inside has a non-slip material to help prevent sliding out of position during the transfer, whilst the outside has a low friction side to aid placement and removal.



SmartMove 2Handles



SmartMove 4Handles

*Etac SmartMove has a non-slip material to help prevent sliding out of position.*

### Intended Use

- Support stand up / sit down
- Support moving higher up in bed
- Moving back / forward in chair / wheelchair
- Support transfer between bed, wheelchair, chair, etc



Description	Width	Length	Item no.
SmartMove 2Handles	21 cm 8.3"	60 cm 23.6"	IM425
SmartMove 4Handles	21 cm 8.3"	60 cm 23.6"	IM426

### Specifications

Outside: Nylon with polyurethan coating  
 Inside: Polyester  
 Handles: Polyester  
 Underside: Polyester



### Washing Instructions

Do not use fabric softener or bleach.

