

Etac SupportBelt

Support belt with handles

Etac SupportBelt is used to offer support and assistance during mobilization and transfers. The diagonally placed handles allow the caregiver to maintain a more ergonomic grip.

SupportBelt is a useful tool when helping a person to stand up or sit down, to support walking or to assist during transfer between bed, wheelchair, chair, etc. SupportBelt has non-slip material on the inside and an adjustable belt buckle.



Etac SupportBelt has diagonally placed handles that allows the caregiver to maintain a more ergonomic grip.

Intended Use

- Support stand up / sit down
- Support walking training
- Support transfer between bed, wheelchair, chair, etc
- Transfer from wheelchair to toilet

Description

SupportBelt XS, 4 handles, waist 65-105 cm
SupportBelt S, 5 handles, waist 65-105 cm
SupportBelt M, 6 handles, waist 90-130 cm
SupportBelt L, 8 handles, waist 125-150 cm
SupportBelt XL, 11 handles, waist 145-170 cm

Width

9 cm 3.5"
13 cm 5.1"
13 cm 5.1"
13 cm 5.1"
15 cm 5.9"

Length

105 cm 41.3"
105 cm 41.3"
130 cm 51.2"
150 cm 59.1"
170 cm 66.9"

Item no.

IM8010509
IM8010513
IM8013013
IM8015013
IM8017015

Specifications

Outside: Nylon
Inside: Nylon with PU-coating
Padding: Polyetylen
Handles: Polyester strap with polypropylen
Belt Buckle: Acetal, nylon

Colour code - size labelling

	XS
	S
	M
	L
	XL

Washing Instructions

Can be disinfected with LT-CIN isopropylalcohol 72 %).
Do not use fabric softener or bleach.

