

Etac Positioning Wedge

Support and positioning cushions with non-slip material

Etac Positioning Wedge supports the user on their side after turning, when resting or when performing hygiene or wound care, thus replacing having a caregiver in a static holding position.

Positioning Wedges can be combined with Etac 4WayGlide System, SatinSheet System, Multi-Glide or other slide sheets. The cover is a non-slip Polyurethane material that helps to keep the Wedge in place independent of surface material.

The Etac Positioning Wedges have unequally angled sides, turn the wedge around and use the angle that fits best for the situation.



Description	Width	Height	Length	Item no.
Positioning Wedge Small	30 cm 11.8"	20 cm 7.8"	53 cm 20.9"	IM53NS
Positioning Wedge Large	30 cm 11.8"	20 cm 7.8"	80 cm 31.5"	IM53/80NS

Washing Instructions

Can be wiped with disinfectant
Do not use fabric softener or bleach.

Specifications

Soft foam
Nylon with PU-coating, zipper



Intended Use

- Support static positioning
- Resting
- Care and personal hygiene
- Can be used alone or in combination with other positioning tools

