EVOCARE CONTACT:

Sales: 07 3355 8000 Fax: 07 3355 5043

Email: sales@evocare.com.au
Web: www.evocare.com.au

PO Box 1144, Stafford Q 4053 Australia

EV3160 X-Ray Viewers

Call or email for more information

More information on our website: http://www.evocare.com.au/x-ray_viewers.html





Codes & Dimensions:	
EV31600 Single:	377mm x 500mm x 105mm
EV31601 Double:	750mm x 500mm x 105mm
EV31602 Triple:	1125mm x 500mm x 105mm
EV31603 Quadruple:	1500mm x 500mm x 105mm

Flush Mount Models

All of the above surface mount models are available as flush mount. To specify flush mount, add **FM** after the model number. Flush mount models have a 25mm flange around the front face.

Features:

- · Bench and Flush Mount Wall Models
- Excellent light diffusion and uniform illumination, achieved from:
 - circular fluorescent tubes
 - continuous opal Perspex viewing panel
 - · internal white reflector surface
- Grip rollers in stainless steel strip to hold film
- Individually switched bays internally divided
- Easy to clean white powder coated steel cases
- Illuminated panel size:
 - EV31600 single bay length 335mm x height 400mm
 - EV31601 double bay length 710mm x height 400mm
- Viewer Box depth is 105mm

