

GEL PRODUCTS


Capillary closure is at the root of ulcer formation.

The blockage or damage of capillaries prevents the natural maintenance of healthy skin.

In operation, the use of proper positioning devices prevent excessive pressure build-up and capillary closure.

Evenly distribute weight and provide constant support beneath patients.

Gel products are cost effective, latex free, fast & easy cleaning, radiolucent tissue equivalent.



Operating room

DOUGHNUT GEL HEAD PAD



SK-22

Doughnut gel head pad provides excellent protection when the patient is positioned for surgery by supine position.

SK-22 (Adult) 20(D) * 4.5(H) cm

C-DOUGHNUT GEL HEAD PAD



SK-22-1

C-doughnut gel head pad provides clear air passage way when the patient is positioned for surgery by lateral or prone position.

SK-22-1 (Adult) 21(D) * 4.5(H)cm

PRONE GEL HEAD PAD



It provides clear air passage way when the patient is positioned for surgery by prone position.

SK-22-2 27(L)*24(W)*15(H) cm

GEL BODY PAD



SK-23-2

SK-23

SK-23-1

Gel body pad provides complete pressure reduction and shares across entire table surface.

SK-23 160(L) * 50(w) * 1.3 (H)cm

SK-23-1 183(L) * 50(w) * 1.3 (H)cm

SK-23-2 84 (L) * 50(w) * 1.3 (H)cm

GEL CHAIR PAD




Gel chair pad is used underneath prominence area in a variety of applications.

SK-24 42(L)*42(W)*1.3(H)cm

SK-24-2 36(L)*36(W)*1.3(H)cm

SK-24-3 30(L)*30(W)*1.3(H)cm

GEL EMBO PAD



SK-24-1

Gel embo chair pad is used for iliac crest, under heels and elbows. Is used as kidney rest, axillary roll while patient in lateral position.

Gel embo chair pad
 SK-24-1 40(L)*45(W)*1.6 (H)cm

GEL KNEE PAD




Gel knee pad is used for OBGY surgery in lithotomy position. It provides against nerve injury.

SK25 40(L)*24(W)*1(H)cm

Operating room

GEL LATERAL PAD/GEL



SK-26

GEL THYROID PAD



SK-26-2

Gel lateral pad is designed to accommodate both the shoulder and provides exceptional pressure relief and nerve protection.

Gel thyroid pad is designed for thyroidectomy. It provides patient with gentle pressure while long surgery with hyperextension position. The tissue friendly feels make it very comfortable.

SK-26	70(L) * 50(W) * 12(H) / 0.8cm
SK-26-2	45(L)*30(w)*7(H) cm

GEL SHOULDER PAD



SK-27



SK-27-2 (hard)



SK-27-3



SK-27-1 (hard)



SK-27-4

Semi-circular design provides maximum pressure protection while patient is in supine position.

sk-27, 27-4 can be used under the knee, the foot against the board.

While patient is in lateral position, sk-27 can be used for the shoulder protection.

In thyroidectomy, sk27-4 is used for neck hyperextension.

while patient is in prone position sk-27-1, sk-27-2 can be used under the chest.

SK-27	47(L) * 10(w) * 5 (H)cm
SK-27-1	47(L) * 15(w) * 8 (H)cm
SK-27-2	47(L) * 15(w) * 12(H)cm
SK-27-3	25(L) * 15(w) * 8(H)cm
SK-27-4	52(L) * 10(w) * 8(H)cm

GEL HEEL PAD



Gel heel pad provides safety and protects heel area body prominence.

SK-28	20(L)*11(W)*7(H)cm
-------	--------------------

GEL ARMBOARD PAD



It is applied for supine position for protection of bony prominence area.

SK-29	55(L)*15(W)*1.3(H)cm
-------	----------------------