



Strong, reliable and powerful

- SWL 155 or 205 kg
- Designed for intensive use in professional environments
- Meets and exceeds all current norms and regulations in industry
- Strong 24V/4,5AH NiMH battery ensures long usage time, even in heavy use environments
 - Shorter charging time to full capacity
 - Prolonged battery life
 - More lifts from a single charge

- Strong, noiseless and fast Linak actuators
- Audio and visual indicators for low battery
- Varied charging options on board lift and off board lift via extra battery
 - Built in and easily accessible docking station for battery
 - Extra battery packs as option
 - Fully automatic switch mode charger

Easy to maintain

- Removable and easy to clean foot tray
- · Visual service indicator
- Prepared for future software upgrades
- Access to statistics via Service Information Console (SIC) for certified technicians
 - Number of lifts lifespan and since last inspection
 - Number of battery warnings
 - Number of errors/faults

V. Guldmann A/S

Graham Bells Vej 21-23A DK-8200 Aarhus N Tel. +45 8741 3100 Fax +45 8741 3131

Fax +45 8741 3131 info@guldmann.com www.guldmann.com









 an easy, comfortable and safe way to raise and support patients with partial or complete control of upper body

One lift for many functions

- · Sit to stand
- Standing balance (with foot tray removed)
- · Transfers from chair to chair
- · Transfers from bed to chair
- · Transfers from bed to toilet

Adjusts easily to comfortably fit most heights and sizes

- · Padded and adjustable knee support
- · Multiple sling attachment points
- Ergonomically profiled push handle for all caregiver sizes
- Multiple sling sizes to create optimum comfort and support (option)
- Add-on hip sling to assist with hips and buttocks raising (option)

Easy and safe to use and understand

- Electric up, down, and leg width adjustment
- Dual operations options hand control or control box
- Rear foot push bar to increase maneuverability
- Manually and electrically operated emergency lowering
- · Easily accessible emergency stop
- User guide permanently attached to control box – 'How to' review always available