

EVO CARE CONTACT:

Sales: 07 3355 8000

Fax: 07 3355 5043

Email: sales@evocare.com.au

PO Box 1144, Stafford


Q 4053 Australia

LMT5 - TENS & EMS

For the treatment of facial muscles in stroke victims

Call or email for more information

Website Link: <http://www.evocare.com.au/tens-lmt5.php>



Side View



LMT5 - Digital TENS (*Transcutaneous electrical nerve stimulation*) and EMS (*Electric Muscle Stimulation*) with Soft Bag

Technical Specifications :

Digital Buttons:	Control EMS and TENS modes with large LCD display.
Dual Channel:	Isolated between channels.
Modes:	3 TENS and 3 EMS
Pulse Rate:	Adjustable. 1Hz - 150Hz
Pulse Width:	Adjustable. 30 – 300 microseconds
Pulse Intensity	Adjustable 0-100mA peak into 500ohm load each channel
Timer:	5 to 95 minutes and Continuous.
Wave Form:	Asymmetrical Bi-phasic square pulse
Ramp Time:	1—20 seconds.
Cycle On Time:	On: 1—30 seconds, Off: 1—60 seconds.
Battery:	3 x AAA batteries (4.5V)
Warranty:	1 year

Package contains:

- LMT5+EMS TENS unit with Belt Clip
- Carry Case
- LMTEL3200 3.2cm diameter electrodes (pack 4) self adhesive and pre-gelled
- 5cm square electrodes (pack 4) self adhesive and pre-gelled
- Set of Leads
- 3 x AAA Batteries
- Instruction Booklet



Top View

Available Replacement Electrodes:

LMTEL5090 Electrode (Pack 4) 50 x 90mm square self adhesive and pre-gelled
LMTEL Electrode (Pack 4) 50 x 50mm square self adhesive and pre-gelled



Our catalogue is our comprehensive and up-to-date website: www.evocare.com.au

Specialist Equipment Sales & Service