



# Anaesthetic Apparatus DRYING CABINETS



 **Evocare**<sup>®</sup>  
www.evocare.com.au

Phone: 07 33558000 Fax: 07 33555043



## FLEXIBILITY

A cabinet which can handle large loads with any combination of tubes and parts.

## EASY INSTALLATION

Connects to standard hospital power outlets; a vent connection removes heated air.

## COMPACT

Centre opening doors protrude only 330mm when fully open. Pass through wall units are available.

## ECONOMY

Double wall stainless steel construction with an insulated cavity, combined with double glazed doors. Heating is discontinued when doors are open.

## COMPLIANCE

Meets both AS2514 (Drying cabinets for Medical Equipment) and AS2774 (Drying cabinets for Respiratory Equipment). Electrical Safety CS071150V.

## APPLICATION

The Malmel DC series of Drying Cabinets are designed to dry respiratory tubing, face masks, surgical instruments and ancillary equipment.

## SAFETY

Dual over temperature cut-out. Set point temperatures are locked to prevent unauthorised change. Cabinet temperature is pre-set at 70°C and governed not to exceed this temperature.

## VERSATILITY

A modular hanging system with shelving and trays to suit any requirement.

## EFFICIENCY

A centrifugal fan drives heated air downwards to provide rapid internal drying of suspended corrugated tubing and circulate warmth throughout the chamber.

## HYGIENE

A washable filter removes fluff, lint, dust etc.



 **Evocare**<sup>®</sup>  
www.evocare.com.au

Phone:07 33558000 Fax:07 33555043

\* ACCESSORIES AND PRODUCTS SHOWN IN ALL PHOTOS ARE NOT INCLUDED WITH UNITS.

# MALMET

## INSTRUMENT AND ANAESTHETIC APPARATUS DRYING CABINET

- Handles large loads with any combination of instruments, tubes, masks and parts
- Designed to dry surgical instruments, respiratory tubing and face masks
- Connects to standard hospital power outlets
- Compact design
- Double wall stainless steel construction
- Double glazed doors
- Dual over temperature cut-out and cabinet preset temperature of 70°C
- Available also as "pass through" cabinet
- ARTG Certification 188763 Class I

<b>UNIT DIMENSIONS:</b>	All front loading cabinets	750mm (W) x 640mm (D) x 1750mm (H)
	All "pass through" cabinets	750mm (W) x 690mm (D) x 1750mm (H)
<b>VOLTAGE:</b>	240V	
<b>CURRENT:</b>	15 Amps	
<b>PHASE:</b>	Single	
<b>HERTZ:</b>	50	
<b>SHIPPING WEIGHT:</b>	185 kg	
<b>VOLUME:</b>	475 litres	
<b>SET TEMPERATURE:</b>	70°C	
<b>POWER CONSUMPTION:</b>	Less than 78kW/24hrs	
<b>ACCESSORIES</b>	Anaesthetic Bag Holder (Product Code 87-5038) Tube Holder Rack (Product Code 87-5104) Full Width Wire Shelf (Product Code 87-5139) Half Width Wire Shelf (Product Code 87-5215) Door Lock Bolt Assembly; fitted in factory (Product Code 100-5506)	



MODEL:DC12



### FRONT LOADING DRYING CABINET

Product Code	Model No	Description
8-A00-87-030	DC11	MALMET DRYING CABINET – FRONT LOAD (6 HALF WIDTH WIRE SHELVES; 2 RACKS OF 18 x 37MM DIA TUBE HOLDERS)
8-A00-87-035	DC11T	MALMET DRYING CABINET – FRONT LOAD (6 HALF WIDTH WIRE SHELVES; 2 RACKS OF 18 x 37MM DIA TUBE HOLDERS; WITH TIMER)
8-A00-87-040	DC12	MALMET DRYING CABINET – FRONT LOAD (6 FULL WIDTH WIRE SHELVES; NO TUBE HOLDERS)
8-A00-87-045	DC12T	MALMET DRYING CABINET – FRONT LOAD (6 FULL WIDTH WIRE SHELVES; NO TUBE HOLDERS; WITH TIMER)
8-A00-87-051	DC13	MALMET DRYING CABINET – FRONT LOAD (FULL HANGING; 5 RACKS OF 18 x 37MM DIA TUBE HOLDERS)
8-A00-87-055	DC13T	MALMET DRYING CABINET – FRONT LOAD (FULL HANGING; 5 RACKS OF 18 x 37MM DIA TUBE HOLDERS; WITH TIMER)

### "PASS THROUGH" DRYING CABINET

Product Code	Model No	Description
8-A00-87-060	DC21	MALMET DRYING CABINET – FRONT LOAD & REAR UNLOAD (6 HALF WIDTH WIRE SHELVES; 2 RACKS OF 18 x 37MM DIA TUBE HOLDERS)
8-A00-87-065	DC21T	MALMET DRYING CABINET – FRONT LOAD & REAR UNLOAD (6 HALF WIDTH WIRE SHELVES; 2 RACKS OF 18 x 37MM DIA TUBE HOLDERS; WITH TIMER)
8-A00-87-070	DC22	MALMET DRYING CABINET – FRONT LOAD & REAR UNLOAD (6 FULL WIDTH WIRE SHELVES; NO TUBE HOLDERS)
8-A00-87-075	DC22T	MALMET DRYING CABINET – FRONT LOAD & REAR UNLOAD (6 FULL WIDTH WIRE SHELVES; NO TUBE HOLDERS; WITH TIMER)
8-A00-87-081	DC23	MALMET DRYING CABINET – FRONT LOAD & REAR UNLOAD (FULL HANGING; 5 RACKS OF 18 x 37MM DIA TUBE HOLDERS)
8-A00-87-085	DC23T	MALMET DRYING CABINET – FRONT LOAD & REAR UNLOAD (FULL HANGING; 5 RACKS OF 18 x 37MM DIA TUBE HOLDERS; WITH TIMER)

<b>Accessories / Spare Parts</b>	
<b>PRODUCT CODE</b>	<b>DESCRIPTION</b>
8-M00-87-033	Instrument Tray, 362 x 271 x 22mm
8-M00-87-034	Air Diffuser Blank
8-M00-87-035	Air Diffuser (with holes)
87-5038	Anaesthetic Bag Holder, Double Sided 6
87-5215	Part width shelf DC11/DC21 models. Size 520d x 340w
87-5139	Part width shelf DC12/DC22 models. Size 520d x 600w
8-M00-87-095	10 x 50mm Diameter tube holder
8-M00-87-100	16 x 37mm Diameter tube holder (models prior to Feb 2002)
87-5104	18 x 37mm Diameter tube holder
100-5506	Door Lock Bolt Assembly (must be factory fitted)

Queensland Distributor

# EVOCARE AUSTRALIA PTY LIMITED

A.B.N. 98 078 566 604

Trading as **EVOCARE** and **L&M EQUIPMENT**

P.O. Box 1144, Stafford Qld. 4053

Ph: 07 3355 8000 Fax: 07 3355 5043

Website: <http://www.evocare.com.au>

Email: [sales@evocare.com.au](mailto:sales@evocare.com.au)

[workshop@evocare.com.au](mailto:workshop@evocare.com.au)

[warehouse@evocare.com.au](mailto:warehouse@evocare.com.au)

[accounts@evocare.com.au](mailto:accounts@evocare.com.au)



Global-Mark.com.au®

## Product Range

Tall Warming Cabinet  
**FWTFC**



Malmet Instrument and Anaesthetic Apparatus Drying Cabinet  
**DC12**



Malmet Fluid and Blanket Warming Cabinet  
**FBWFC**



Malmet Blanket Warming Cabinet - with Stand  
**BWS**



Malmet Utensil/Bowl Washer Disinfector, Top Loading  
**WDT2**



Malmet Bedpan/Urinal Bottle Washer Disinfector  
**ES-D**



Malmet Bedpan/Urinal Bottle and Utensil/Bowl Washer Disinfector  
**WDS**



Malmet Rack (for Bedpan & Urinal Bottles with Wall Plate & Drip Tray)  
**BR6-PR6 + WT**

## National Head Office & Customer Service

### For Correspondence:

9-11 McKay, Leeton, NSW 2705 Australia

### General Sales, Service, Export Enquiries & Accounts

Ph. +61 2 6953 7677 Email: [info@malmet.com.au](mailto:info@malmet.com.au)

### Warning:

This brochure contains general information only. If you require accurate dimensions or any other details which may affect your installation, please contact Malmet. E.&O.E. In the interest of continued product improvement, Malmet reserves the right to alter design and product specifications without notice. E.&O.E. For further information on warranty terms and conditions, please refer to Malmet's Warranty Statement.

MEDICAL RECEPTACLES SHOWN IN ALL PHOTOS ARE NOT INCLUDED WITH UNITS.

PBL005 0119

