

MALMET

BEDPAN/URINAL BOTTLE AND UTENSIL/BOWL WASHER DISINFECTOR

The WDS Series combines the performance, capability and proven reliability of Malmets ES-D Series Bedpan/Urinal Bottle Washer Disinfector and WDT Series Utensil/Bowl Washer Disinfector. This offers great space and cost savings for your facility. Being locally made, the availability of spare parts and servicing is a simple and quick process. Malmets WDS has been micro-biologically tested to verify that the machine effectively inactivates organisms (bacteria and viruses) as required under the relevant Australian and ISO standards.

- Only combination washer disinfector designed and manufactured in Australia
- Large capacity: 2 bedpans and 4 urinal bottles, or up to 4 bowls, or up to 6 kidney dishes (with optional rack), various utensils etc
- Infection control: Microbiologically tested by an independent NATA certified laboratory
- Hands-free operation to minimise infection control risk
- Sensors to minimise risk of incorrect program being selected
- Simple operation: 3 programs, with clear graphics (bedpans/urinal bottles, urinal bottles only, or utensils/bowls)
- All programs fully disinfect
- Detergent wash cycle
- Steam disinfection with inbuilt boiler
- Disinfection to higher A₀600 level, as recommended by AS4187:2014/ISO15883
- No venting required
- Available in 415V 20A or 240V 20A
- Compliant to relevant electrical safety and Watermark standards
- ARTG Certification 232450 to industry highest Class IIb level
- The door obstruction feature is factory set; causing the door to re-open when an obstruction is detected. This is designed to protect the integrity of the unit and the articles within.



MODEL: WDS



UNIT DIMENSIONS:	605mm (W) x 625mm (D) x 1685mm (H)		
VOLTAGE:	240V or 415V		
PHASE:	1 phase or 3 phase		
HERTZ:	50		
SHIPPING WEIGHT:	165 kg		
CAPACITY:	2 small or 2 bariatric slipper pans and 4 urinal bottles OR 2 large bowls OR 3 medium bowls OR 4 small bowls Any combination of utensils (tooth mugs, kidney dishes etc)		
DETERGENT BOTTLE CAPACITY:	5 Litre (Malmets Product Code 81-9071SP) Detergent ARTG Certification 190013		
DETERGENT DISPENSED PER APPLICATION:	30 ml		
DETERGENT SPECIFICATIONS / GRADE:	Poisons Schedule: 6 Dangerous Goods Class: 8 Hazchem Code: 2R		
1 PHASE MODELS WASH SETTINGS:	Bedpan/Urinal bottle	Urinal bottle only	Bowl/Utensils
AVERAGE CYCLE TIME:	15.15 minutes	14.25 minutes	14.25 minutes
POWER CONSUMPTION:	0.62 kW/hr	0.572 kW/hr	0.572 kW/hr
WATER CONSUMPTION:	35.6 litres (cold water) 21.6 litres (hot water)	24.4 litres (cold water) 20.7 litres (hot water)	24.4 litres (cold water) 20.7 litres (hot water)
3 PHASE MODELS WASH SETTINGS:	Bedpan/Urinal bottle	Urinal bottle only	Bowl/Utensils
AVERAGE CYCLE TIME:	14.15 minutes	13.15 minutes	13.15 minutes
POWER CONSUMPTION:	0.60 kW/hr	0.574 kW/hr	0.574 kW/hr
WATER CONSUMPTION:	35.6 litres (cold water) 21.6 litres (hot water)	24.4 litres (cold water) 20.7 litres (hot water)	24.4 litres (cold water) 20.7 litres (hot water)
ACCESSORIES:	Kidney Dish Rack (Product Code 92-5033), Bottle Holder for 1 litre jug (Product Code 92-5151), Jug Holder Rack for 3 x 1 litre jugs (Product Code 92-5153), Utensil Basket (Product Code 92-5235)		

BOWL RACK INCLUDED HOLDS- UP TO 4 BOWLS



OPTIONAL: UTENSIL BASKET FOR SMALLER ITEMS

AVAILABLE MODELS	Product Code	Model No	Description
	92-1000	WDS1S	MALMET BEDPAN/URINAL BOTTLE AND UTENSIL/BOWL WASHER DISINFECTOR (S TRAP, 1 PHASE, 20 AMP)
	92-1001	WDS1P	MALMET BEDPAN/URINAL BOTTLE AND UTENSIL/BOWL WASHER DISINFECTOR (P TRAP, 1 PHASE, 20 AMP)
	92-1010	WDS3S	MALMET BEDPAN/URINAL BOTTLE AND UTENSIL/BOWL WASHER DISINFECTOR (S TRAP, 3 PHASE, 20 AMP)
	92-1011	WDS3P	MALMET BEDPAN/URINAL BOTTLE AND UTENSIL/BOWL WASHER DISINFECTOR (P TRAP, 3 PHASE, 20 AMP)