




seca

seca 201










♂	< 102 cm*	≥ 102 cm*
♀	< 88 cm*	≥ 88 cm*

*Obesity: preventing and managing the global epidemic. Report of a WHO Consultation. WHO technical report series 894 (1999: Geneva Switzerland).



• Technical Data

	205 cm
	1 mm
	± 5 mm
	$\varnothing = 70$ mm $h = 22$ mm

	~ 55 g
	10°C – 40 °C
	-25°C – 65°C

• Cleaning

Clean the measuring device as required using a domestic cleaning agent or commercially-available disinfectant. Follow the manufacturer's instructions.

Guarantee

We guarantee defects attributable to poor materials or workmanship for 24 months from the date of purchase.

The guarantee is invalidated if the device is opened by unauthorised parties.

Important: guarantee valid only in conjunction with proof of purchase.

Declaration of Conformity 0123

Mechanical measuring device for general measurement of the circumference of children and adults and for measuring the length of extremities.

Model: seca 201

Batch no.: see grip of tape measure

The measuring device complies with the fundamental requirements of Directive 93/42/EEC pertaining to medical products.

EVOCARE AUSTRALIA PTY LIMITED

A.B.N. 98 078 566 604

Trading as **EVOCARE** and **L&M EQUIPMENT**

P.O. Box 1144, Stafford Qld. 4053

Ph: 07 3355 8000 Fax: 07 3355 5043

Website: <http://www.evocare.com.au>

Email: sales@evocare.com.au

workshop@evocare.com.au

warehouse@evocare.com.au

accounts@evocare.com.au

