

October 2020

Etac Turner Pro

Transfer platform - user turning with standing support

A sit to stand aid with functional design

Etac Turner Pro is a turn aid with a functional design that offers safe patient turning with standing support. This turn aid enables a good posture and minimises the risk of injury for the carer when transferring, moreover, it encourages natural standing and sitting movement for the user. The only adjustment necessary is the height of the leg support and the handle. Thereafter, the transfer can be performed safely whilst the caregiver and user can maintain eye contact and communication.

The orange handle provides a strong, high contrasting colour which can help users with impaired vision or dementia. The small footplate makes it the perfect choice of transfer aid when transferring in tight spaces. Etac Turner Pro is easy to transport due to its low weight of only 7.5 kg, furthermore it is compact and easy to grip.



Description	Dimensions	Item no.
Etac Turner Pro	800 - 1200 mm	16090106
Material		
Cast aluminium, powder co	ated steel, chrome 3	
Weight		
Total weight	165 lbs	
Maximum user weight (SWL	_) 440 lbs	
Dimensions		
Lower leg support height	150 mm	
Lower leg support width	380 mm	
Facts		
	max max	
Alk. pH 70% 5-9 3 min	J 200 440	CE

Features

- The sit to stand aid is easy to transport due to a low weight of only 7.5 kg. It is compact, evenly balanced and easy to grip.
- The angled leg support is soft and comfortable and provides lateral support. Easy to adjust to whatever height is needed.
- The orange handle provides a strong, high contrasting colour which can help users with impaired vision or dementia.
- The cast aluminium base plate with integrated anti slip function is easy to turn. Together with the wheel, Turner Pro works well on different flooring materials including carpets.



Etac AB

info@etac.se www.etac.com



"Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances"

For the latest news and continuously updated product information – visit: www.etac.com

Queensland Distributor

EVOCARE AUSTRALIA PTY LIMITED A.B.N. 98 078 566 604



Trading as EVOCARE and L&M EQUIPMENT P.O. Box 1144, Stafford Qld. 4053 Ph: 07 3355 8000 Fax: 07 3355 5043 Website: http://www.evocare.com.au Email: sales@evocare.com.au workshop@evocare.com.au warehouse@evocare.com.au accounts@evocare.com.au



Etac AB Kista Science Tower SE-164 51 Kista, Sweden Tel +46 371 - 58 73 30 Fax +46 371 - 58 73 90 info@etac.se www.etac.com