carilex⁶

Dual Flex

Dynamic Hybrid Mattress Solution

When Comfort is a Priority!

The Dual Flex is a dynamic replacement mattress system, combining both high performance static foam, and active alternating air cell technology for exceptional comfort and healing. Dual Flex provides a very flexible step-up or step-down mattress system to cater for higher or lower risk patients all in one, cost effective pressure relieving mattress solution.





Whisper Quiet



Alternating Pressure Relief



CLP Therapy Pressure



Low Pressure Alarm



Power Failure Alarm



Patient Transfer Mode



GENERAL FEATURES

| 1 Power/Mute | Turns pump on and mutes alarm |
|---------------------------------|--|
| 2 LED - Power Failure Alarm | Green LED lights up when power is on Flashes orange when power fails Audible and visual power failure alarms to alert caregivers in case of abnormal power outage |
| 3 Comfort Levels | Auto-set, with 5 manual harder and softer settings |
| 4 Pressure Monitoring Indicator | Monitors mattress internal pressure 24 hours to achieve optimal internal pressure and maximum pressure relief |
| 5 Low Pressure Indicator | Flashing orange LED indicates insufficient internal pressure |
| Dynamic Mode | 2:1 alternating therapy, 10 minute cycles |
| CLP Therapy | All air cells maintain constant low pressure support |
| Ultra-Quiet Pump | Less than 28db, for a peacefull sleeping environment |
| Auto-Set | No need to input weight, with ability to manually adjust for patient comfort |
| Transport Mode | Connect the two quick couplings together to retain air in the mattress, can stay inflated for up to 15 hours |
| | |

TECHNICAL DETAILS

Power Outage Support

| Risk/Category | Very High |
|-----------------------------|--|
| Waterlow Score | 20+ |
| Pressure Sore Grading | 4 |
| Safe Working Load | Up to 250kg |
| Compatible Support Surfaces | Dual Flex Mattress Replacement – Single and King Single |
| Travel/Weight | Pump weighs 2.5kg |
| Clinical Settings | Suitable for Hospital, Home Care and Nursing Home environments |
| Warranty | 2 year |

Can be used with or without pump

Distributor

EVOCARE AUSTRALIA PTY LIMITED

A.B.N. 98 078 566 604
Trading as **EVOCARE** and **L&M EQUIPMENT**P.O. Box 1144, Stafford Qld. 4053
Ph: 07 3355 8000 Fax: 07 3355 5043
Website: http://www.evocare.com.au



P.O. Box 1144, Stafford Qld, 4053
Ph: 07 3355 8000 Fax: 07 3355 5043
Website: http://www.avocare.com.au
Email: sales@evocare.com.au
warehouse@evocare.com.au
warehouse@evocare.com.au
accounts@evocare.com.au




carilex®

Dual Flex

Dynamic Hybrid Mattress Solution

When Comfort is a Priority!





Mattress Replacement



Multiple Mattress Sizes



Easy Cleaning



GELTEX® Medical Grade



Cable Management



Made in UK



Alternating Air



Zones

DUAL FLEX Hybrid Mattress Replacement

Dual Flex's high performance foam layer features a zonal area design, including GELTEX inside technology, to provide effective redistribution of pressure on the sacrum, head, shoulder blades and heel areas. Alternating therapy can be quickly and easily applied by connecting the ultra-quiet pump, whilst still providing a peaceful sleeping environment.

| GENERAL FEATURES | |
|--|--|
| GELTEX * Inside Technology - Foam Topper | Foam is clinically proven, breathable, medical grade and temperature resistant |
| Zonal Area Design | Sacrum, head, shoulder blades and heel areas |
| Low Interface Pressures | For improved blood circulation and wound healing |
| Power Outage Support | Can be used with or without a pump |
| Channels for Hose and Cell Connections Built into Foam Base | Avoids pressure points on the surface |
| Hygiene Control | Durable, vapour permeable and water resistant 2-way stretch coverlet to protect the skin from shearing |
| Nylon Cells | TPU laminations for longevity |
| Flexible Contour Design | For profiling beds |
| Hidden Zip | For improved infection control |
| Buckle and Straps | For securing it to the bed frame |
| TECHNICAL DETAILS | |
| Risk/Category | Very High |
| Waterlow Score | 20+ |
| Pressure Sore Grading | 4 |
| Safe Working Load | Up to 250kg |
| Air Cell Configuration | 19 x 3" cells |
| Cycle Times | 10 minutes, 2:1 |
| Available Size/Dimensions | Single 90cm x 200cm x 18cm, King Single 107cm x 200cm x 18cm |
| Securing Mattress to Bed | Buckle Straps |
| Mattress Weight | 8kg, with handles for ease of carrying |
| Clinical Settings | Suitable for Hospital, Home Care and Nursing Home environments |
| Compatible Pumps | Dual Flex |
| Warranty | 2 year |

