



GH3+ Scale module (accessories)

The scale module is very useful anywhere the precise weighing of an individual is needed to monitor his or her changes in weight.

Scale module

- The scale module can be used for weighing patients whose weight needs to be monitored.
- Medical scales can and may only be calibrated by an accredited individual/company, and Guldmann's integrated scales are therefore not medical scales.
- Guldmann's scale module can be calibrated by a technician approved by Guldmann (local technical staff).
- The accuracy of Guldmann's integrated scale is the same as medically certified scales.
- The GH3+ solution complies with the EN/ISO 10535:2007 standard.

Specifications of the scale module

Capacity 0 - SWL
 Accuracy +/- 0.1 % at max load
 Display resolution (d) 0.1 kg (0.22 lbs)
 Accuracy of repetition < 0.1 kg at 0-200 kg (0-440 lbs)
 < 0.2 kg at 0-350 kg (0-770 lbs)
 Minimum weight 5 kg (11 lbs)
 Display type LCD in Hand Control