



Pleto


The mobile overbed table

Thanks to its excellent manoeuvrability, the Pleto server can be used anywhere in everyday care situations. Usually, however, it is used as a table. It can be pushed right up to beds, chairs or recliner chairs, as required. Its convenient C frame guarantees high strength and stability.

With its gas spring support, the server can be adjusted from a height of 69 cm to 104 cm. This makes it possible to adapt the server flexibly to various situations, such as for mealtimes or as an ergonomic bookrest for reading. The table top can be tilted, but a small section of the table always remains horizontal, creating a practical surface on which objects can be placed. This allows patients to put down their drink while reading, for example.

The server height is adjusted automatically when the bed height is adjusted, since upward adjustments are possible without releasing a handle. This is effective in avoiding collisions.

External dimensions/height adjustment range:



◀ Split and horizontal tilting table top



▲ The height adjustment lever is within easy reach



▲ A strong chassis provides stability