

carilex®

CariChair

Semi-automatic mobile cushion system with battery backup – defying the constraint of time and place for uninterrupted, around-the-clock care. Whenever, wherever.





Semi-automatic mobile cushion system with battery backup – defying the constraint of time and place for uninterrupted, around-the-clock care. Whenever, wherever.

Giving you the care you need, whenever, wherever. The Carilex CariChair is the advanced mobile cushion system that offers uninterrupted, active therapy for patients on-the-go. With CariChair, the patients no longer need to be confined to the bed for pressure ulcers treatments. An exceptionally progressive accessory to the extensive Carilex wound care solution family, CariChair is designed to be compatible with a variety of mobile devices including geriatric chair, recliner and wheelchair. Packed with the leading-edge IPS technology, CariChair responds to patient movements and automatically adjust internal cushion pressure. The Semi-Automatic Power Unit does not need to input patient's weight. Thanks to years of rigorous research and countless examinations in what is most comfortable for patients, CariChair offers multiple comfort levels that delivers optimal internal pressure according to each patient's sensation. The rechargeable battery backup keeps the system operating up to 8 hours, allowing more flexibility in its usage. Carilex took all the complications in wound care systems and turn it into carefully calculated technology that is fuss-free for you. The most advanced technology packed in one small device, CariChair is the perfect mobile system for total and complete wound care therapy.

The CariChair Power Unit with IPS Technology.

The IPS Technology, Intelligent Pressure Sensing Technology, responds to patient movements on cushion by automatically adjusting internal cushion pressure, which allows CariChair to regulate airflow and continuously provide therapy for our patients. Carilex is the first brand in the industry to master IPS. This gives us more experience and vast knowledge in the technology that is years ahead of our competitors.

1. Multiple compatibility

Compatible with multiple cushions for various types of chairs.

2. Intelligent pressure sensing

Responds to patient movements on mattress by automatically adjusting internal mattress pressure.

3. Comfort setting

Provides optimal internal pressure according to each patient's sensation. Includes options for comfort level adjustment, no weight input needed.

4. Battery backup upto 8 hours

Automatic power detection sensor activates battery backup, in the even of main power failure and allows treatment to continue uninterrupted during patient transportation and examination.

5. Dynamic / Static therapy modes

Provide both dynamic and static functions to meet individual patient's need.

6. Multiple cycle times

Various cycle time selection for alternating pressure therapy.

The CariChair Cushion – for geriatric chair and recliner

Providing uninterrupted dynamic / static therapy, CariChair fits seamlessly with recliner chairs while offering pressure relief and redistribution for patients with sensitive tissues. Quick release coupling for fast patient transferring and provides an easy connection between power unit and mattress. It supports high weight capacity. CariChair Cushion is the perfect accessory for patients who need to change from bed to chair or vice-versa.

1. Extra support at the sacrum area

Specially designed deep air cells, for excellent immersion and pressure relief.

2. Whole body therapy


Head-to-toe protection for patients even in sitting position.

3. Seamless transition from bed to chair

Continuous, uninterrupted care whenever, wherever.

4. Quick release coupling

Provides an easy connection between power unit and mattress.



Mobicare Plus



Optional Cushion for Wheelchair Paired with Mobicare Plus

1. Extra support at the sacrum area

Specially designed deep air cells, for excellent immersion and pressure relief.

2. Durable material

High quality material intended for long-term usage.

3. Quick release coupling

Provides an easy connection between power unit and mattress.

4. Better immersion

10 3.5" deep air cell.

**Optional cushion for wheelchair
Paired with Mobicare**

1. T-Shape foam base

Specifically designed foam seating for patients to provide better ergonomics.

2. Quick release coupling

Provides an easy connection between power unit and mattress.



Mobicare

Other Features



Battery level LED



Digital LED screen



Whisper Quiet



Panel lock



Alarm Mute



Washable, waterproof cushion cover



Visual and audible indicator



Carry bag



Anti-Skid Security

Technical Data

- Power Unit Model : _____ REF : SR321
- Power Input : _____ 100 - 240Vac 50 60Hz 1.5A
- Power Consumption : _____ 27W Adaptor + Pump; 9W Pump only
- Weight : _____ Power Unit : 0.92kg / 2.03lb
Cushion Set : CariChair cushion : 2.2kg / 4.85lb
MobiCare Plus / MobiCare cushion : 1.5kg / 3.31lb
- Dimension : _____ Power Unit : 24.5 x 15.5 x 6.1cm
CariChair Cushion : 70" x 20" x 5.5" (180 x 50 x 14 cm)
Mobi Care Plus Cushion : 18" x 18" x 3.5" (45 x 45 x 9 cm)
Mobi Care Cushion : 14" x 18" x 3.5" (36 x 45 x 9 cm)
18" x 18" x 3.5" (45 x 45 x 9 cm)
24" x 18" x 3.5" (60 x 45 x 9 cm)
- Agency Approval : _____ SGS-Q Mark UL60101-1
- Operating Condition : _____ Temperature : 5°C to 40°C
Humidity : 15% to 60%
- Transportation & Storage Condition : _____ Temperature : -25°C to +70°C
Humidity : 0% to 93%
- Weight Capacity : _____ 98kg
- This Power Unit is not AP/APG protected

