

# Etac BedString

## String with handles

Etac BedString with handles is intended for the active user by providing support when moving to a sitting position in bed. In the sitting position, the user may use it to lean themselves forward to allow easier application of a sling, positioning of pillows or for dressing.

BedString has 8 handles, alternating right and left side of the string to enable the user to find a good grip.

*Etac BedString with 8 handles on alternating sides of the string.*



### Intended Use

- Support sitting up in bed
- Transfer to bed
- Transfer out of bed



| Description            | Width     | Length       | Item no. |
|------------------------|-----------|--------------|----------|
| BedString with handles | 7.5 cm 3" | 170 cm 66.9" | IM429 ✓  |

### Specifications

Polyester



### Washing Instructions

Do not use fabric softener or bleach.

