


Ravello

Clinic | Care | **Homecare** | Assist



The care chair
for extra mobility

 **Evocare**

www.evocare.com.au
Ph: 07 3355 8000 Fax: 07 3355 5043
Email: sales@evocare.com.au

Model Ravello-Curo "S"

Backrest can be lowered to a 43° reclining position.

Ergonomically shaped seat

The armrests can be lowered to the seat height which makes it easier to stand up and sit down

Wide armrests for a good forearm position

Lever: The reclining position can also be adjusted by the occupant

Rear castors can be locked in pairs

Length-adjustable, fixed-position footrest that tilts to follow backrest adjustments

Safe working load: 200 kg (static), 135 kg (dynamic)

Ravello-Sano and Ravello-Curo - The care chair for extra mobility

Advantages

Safe, comfortable and practical

These are the outstanding features of our Ravello. Advantages which benefit occupants and care staff alike. The upholstered, ergonomically shaped seat and footrests offer maximum sitting and reclining comfort. The front disc wheels make the chair considerably easier to manoeuvre, particularly when the chair is occupied. The Ravello is straightforward to use and genuinely easy to care for with its wipeable seat covers. It is therefore ideal as a care chair.

The basic model for everyday care

The easy-to-operate backrest allows the ergonomically shaped seat to incline to a 43° reclining position, which enhances the comfort of the occupant.

With or without "S"?

The Ravello has 2 different sizes of castor. The rear castors serve as steering castors that can also be locked in position. Of course, the Ravello is also available as an S model. This model is particularly easy to manoeuvre due to its 4 steerable castors (large image).

Design details



Ravello-Sano: The footrest can be fully retracted



Ravello-Curo: Length-adjustable, fixed position footrest that tilts to follow backrest adjustments



Armrest height adjustment in 6 stages



Ravello-Sano and Ravello-Curo are both available in an S version for easier manoeuvrability

Ravello-Sano and Ravello-Curo care chairs

Features

- Backrest: Height: 810 mm, 43° inclinable
- Footrest: 260 x 240 mm, steel tube reinforced plastic, floor clearance 80 mm, non-tipping; anti-slip surface
- Castors: Front castors: Ø 200 mm, fixed
Rear castors: Ø 150 mm, swivel type, can be locked in pairs
- S-version castors: Front castors: Ø 150 mm, swivel type, can be locked in pairs
Rear castors: Ø 150 mm, swivel type, can be locked in pairs
- Seat: PU foam with wooden core, density 50 kg/m³, thickness 80 mm
- Backrest: PU foam with steel core, anatomically formed, headrest with adjustable neck cushion
- Cover material: Imitation leather, seat sewn and bonded, fully water-repellent



We reserve the right to change technical specifications without prior notice.
Last updated: 01/2016, Version 3, Item no. 600201.

Accessories

(available at extra cost)



Side supports



Special headrest with lateral support



Urine bag holder



Attachable table



Infusion stand

Order numbers and external dimensions

30.6439.1 **Ravello-Curo** care chair, overall width: 680 mm, overall height: 1260 mm, overall depth: 1030 mm, seat: 520 x 520 mm

30.6439.2 **Ravello-Sano** care chair, overall width: 680 mm, overall height: 1260 mm, overall depth: 1030 mm, seat: 520 x 520 mm



BURMEIER

The Stieglmeyer-Group

Burmeier GmbH & Co. KG

Pivitsheider Straße 270, 32791 Lage, Germany



www.evocare.com.au

Ph: 07 3355 8000 Fax: 07 3355 5043

Email: sales@evocare.com.au